

ADVANCING HEALTHCARE TECHNOLOGY



PRODUCT CATEGORY:
HOSPITAL FURNITURE

PRODUCT NAME:
HU-I-U0088

PRODUCT DESCRIPTION:
DELIVERY BED - ELECTRICAL



Area of Use: Delivery Rooms / Maternity Wards / Examination Rooms

Intended Use: The bed is intended for obstetrics and gynaecological patients, in all stages of labour, delivery and confinement. These beds have special features both for the comfort and well-being of the patient and for the convenience of healthcare workers.



HU-I-U0088 SPECIFICATIONS

SPECIFICATIONS:	
Motors:	Controlled by 3 x electric motors.
Head/Foot Board Material: (Detachable)	High Density Poly Ethylene head & foot board x 1 set.
Safety/Side Rails:	High Density Poly Ethylene, half length, tuck away side rail x 1 set.
Angle Indicator:	Angle indicator on the side rails.
Bumper:	Bumper wheels on each corner of the bed.
Castors:	4 x castors with central locking system. 150mm diameter
Controllers:	One hand held controller to adjust the movements.
Frame Material:	Epoxy coated mild steel frame.
Minimum Height:	400 mm
Maximum Height:	900 mm
Backrest Angle:	0 - 70°
Trendelenburg / Reverse-Trendelenburg Position:	12° / 12°
Provision for IV Pole:	4 x IV Pole holes provided.
Handles:	CPR handles on each side of the bed.
Optional:	IV Pole
Safe Working Load:	230 Kg

SUPPLIED WITH:	
Waste Bowl:	1 x stainless steel bowl.
Lithotomy Poles:	1 x set Lithotomy Poles
Calf Support:	1 x set fully padded Calf Support
Knee Crutches:	1 x set Knee Crutches
Stirrups:	1 x set Stirrups
Dustproof Cover:	1 x Bottom with dust proof cover.
Mattress:	1 x 3 Section Mattress
Mattress Cover:	1 x Water resistant, antibacterial and breathable mattress cover.