

OUR TRUSTED LEGACY IS YOUR SUCCESS



PRODUCT CATEGORY:
HOSPITAL FURNITURE

PRODUCT NAME:
HU-I-U0002152

PRODUCT DESCRIPTION:
DRESSING TROLLEY



Area of Use: All Areas

Intended Use: Designed to transport and organise medical supplies required for wound dressing procedures within a hospital and clinical environment.

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-I-U0002152 SPECIFICATIONS

TECHNICAL SPECIFICATION:

Dressing trolley with two flat, non-removable shelves.
Heavy-duty carriage mounted on four swivel castors, two of which have brakes.
Transfer bars connect all lower distal portions of the four supports, providing maximum structural strength.
Shelves have a smooth finish, allowing for easy cleaning.
Push-bar handle fitted on one short side.
Both braked castors are mounted on the push-bar handle side.
Top and bottom shelves have guard rails along one length and both widths.
Fitted with protective bumpers at all four corners.

MATERIALS:

High resistance to corrosion (tropical environment).	
Frame and tray:	Austenitic stainless steel 18/10 (304 quality).
Caster frame / bracket:	Steel.
Caster brake:	Total-lock type (wheel and rotational lock).
Caster wheel:	Single wheel, mould-on type, non-hooded (for easy maintenance).
Wheel bearing:	Sealed bearing in the swivel and the wheel.
Swivel:	Ball bearing.
All furniture components contain less than 100 parts per million of the following chemical groups:	
Urea formaldehyde	
Heavy metals including mercury, cadmium, lead, and antimony	
Hexavalent chromium in plated finishes	
Stain and non-stick treatments derived from perfluorinated compounds (PFCs), including perfluorooctanoic acid (PFOA)	
Added antimicrobial treatments.	

DIMENSIONS:

Overall:	90 x 55 x 100 cm (L x W x H)
Frame:	3.2 cm (outside, across), 1.5 mm (thickness)
Shelf thickness:	1.5 mm
Swivel castor wheel:	3.1 x 10 cm (W x diameter) (±5%)
Carrying capacity:	100 kg
Knock-down construction:	Yes

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-I-U0002152 SPECIFICATIONS

SUPPLIED WITH:

- 1 x Complete set of tools required for assembly.
- List of accessories and parts.
- Detailed step-by-step instructions for assembly and safe use, text-free and pictorial-based (i.e. line drawings only).
- Instructions for use, including cleaning and disinfection instructions, in English, French, and Spanish.

PACKAGING, LABELLING, INSTRUCTIONS:

- One (1) unit per box
- Unit weight: 23.3 kg (including packaging)
- Unit volume: 0.06 m³ (including packaging)
- Symbols in accordance with ISO 15223.
- Packaging made of recycled or recyclable materials.