

OUR TRUSTED LEGACY IS YOUR SUCCESS



PRODUCT CATEGORY:
HOSPITAL FURNITURE

PRODUCT NAME:
HU-I-U0002154

PRODUCT DESCRIPTION:
ADJUSTABLE STOOL



Area of Use: All Areas

Intended Use: The adjustable stool is intended to provide seated support for healthcare professionals as an alternative to prolonged standing during clinical procedures, examinations, and patient care activities within hospital and healthcare environments.

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-I-U0002154 SPECIFICATIONS

TECHNICAL SPECIFICATION:

- Movable, height-adjustable stool with a pre-assembled gas spring.
- Heavy-duty carriage mounted on five swivel castors.
- The height-adjustment lever is conveniently positioned beneath the seat.
- The seat is fitted with thick upholstery.
- A safety brake prevents exceeding the maximum height setting.

MATERIALS:

Designed for high resistance to corrosion (tropical environments).

Support and telescopic column:	Chrome-plated
Carriage base:	Single moulded unit, polypropylene or nylon
Castor frame / bracket:	Polypropylene or nylon
Castor wheels:	Twin-wheel, solid type, non-hooded (for easy maintenance), polypropylene or nylon
Wheel bearings:	Polypropylene
Castor fixation into carriage base:	Stainless steel
Upholstery:	High-density polyurethane foam, density 27–33 kg/m ³
Cover:	Flexible plastic cover that is highly tear-resistant, anti-static, flame-retardant, non-absorbing, waterproof, and cleanable with hospital-grade disinfection products.

All furniture components contain less than 100 parts per million (ppm) of the following chemical groups:

Urea formaldehyde.

Heavy metals, including mercury, cadmium, lead, and antimony.

Hexavalent chromium in plated finishes.

Stain- and non-stick treatments derived from Perfluorinated Compounds (PFCs), including Perfluorooctanoic Acid (PFOA).

Added antimicrobial treatments.

DIMENSIONS:

Stool, overall:	60 × 60 × 41–54 cm (L × W × H)
Seat diameter:	35 cm
Upholstery thickness:	4 cm
Carriage base diameter:	60 cm
Support column (with pre-assembled gas spring):	5 cm (outside diameter)
Telescopic upper part:	2.8 cm (outside diameter)
Swivel castor twin-wheel:	4.2 × 5 cm (width × diameter)
Tolerance:	±5%
Carrying capacity:	150 kg
Knock-down construction:	Yes

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-I-U0002154 SPECIFICATIONS

SUPPLIED WITH:

- One (1) complete set of tools required for assembly.
- List of accessories and parts.
- Detailed, step-by-step assembly and safe-use instructions, text-free and pictorial (line drawings only).
- Instructions for use, including cleaning and disinfection instructions, in English, French, and Spanish.

PACKAGING, LABELLING, AND INSTRUCTIONS:

- One (1) unit per box
- Unit weight: 6.7 kg (including packaging)
- Unit volume: 0.06 m³ (including packaging)
- Symbols in accordance with ISO 15223.
- Packaging made from recycled or recyclable materials.