

OUR TRUSTED LEGACY IS YOUR SUCCESS



PRODUCT CATEGORY:  
**HOSPITAL FURNITURE**

PRODUCT NAME:  
**HU-I-U0002156**

PRODUCT DESCRIPTION:  
**MAYO TROLLEY**



**Area of Use:** All Areas

**Intended Use:** The Mayo stand is intended to support and hold surgical instruments, medical devices, and sterile supplies required during surgical or examination procedures. It is designed to provide a stable, adjustable, and easily accessible surface for healthcare professionals, allowing efficient and safe handling of instruments within the sterile field. The stand is portable, enabling convenient positioning close to the patient in operating theatres, procedure rooms, and other clinical environments.

OUR TRUSTED LEGACY IS YOUR SUCCESS



## HU-I-U0002156 SPECIFICATIONS

### TECHNICAL SPECIFICATION:

Movable, height-adjustable instrument table, Mayo type.

Heavy-duty carriage mounted on four swivel castors, of which two are fitted with brakes and two with anti-static wheels.

The support column is positioned side-on the base, facilitating positioning under low-clearance tables (e.g. operating tables).

Both braked castors are mounted on the short side of the base.

A solid manual lever allows adjustment of the telescopic upper section to the required height.

A safety brake prevents exceeding the maximum height setting.

The upper section is fitted with a smooth-finished, removable instrument tray.

### MATERIALS:

Designed for high resistance to corrosion (tropical environments).

|                        |                                                |
|------------------------|------------------------------------------------|
| <b>Frame and tray:</b> | Austenitic stainless steel 18/10 (304 quality) |
|------------------------|------------------------------------------------|

|                                |       |
|--------------------------------|-------|
| <b>Castor frame / bracket:</b> | Steel |
|--------------------------------|-------|

|                      |                                         |
|----------------------|-----------------------------------------|
| <b>Castor brake:</b> | Total-lock type (wheel and swivel lock) |
|----------------------|-----------------------------------------|

|                      |                                                                                        |
|----------------------|----------------------------------------------------------------------------------------|
| <b>Castor wheel:</b> | Single-wheel, mould-on type, non-hooded (for easy maintenance), anti-static (all four) |
|----------------------|----------------------------------------------------------------------------------------|

|                       |                                                 |
|-----------------------|-------------------------------------------------|
| <b>Wheel bearing:</b> | Sealed bearing in both the swivel and the wheel |
|-----------------------|-------------------------------------------------|

|                          |              |
|--------------------------|--------------|
| <b>Swivel mechanism:</b> | Ball-bearing |
|--------------------------|--------------|

All furniture components contain less than 100 parts per million (ppm) of the following chemical groups:

Urea formaldehyde

Heavy metals, including mercury, cadmium, lead, and antimony

Hexavalent chromium in plated finishes

Stain- and non-stick treatments derived from Perfluorinated Compounds (PFCs), including Perfluorooctanoic Acid (PFOA)

Added antimicrobial treatments

### DIMENSIONS:

|                            |           |
|----------------------------|-----------|
| <b>Height, adjustable:</b> | 75–112 cm |
|----------------------------|-----------|

|                    |                              |
|--------------------|------------------------------|
| <b>Upper tray:</b> | 60 × 40 × 2.8 cm (L × W × H) |
|--------------------|------------------------------|

|                        |        |
|------------------------|--------|
| <b>Tray thickness:</b> | 1.5 mm |
|------------------------|--------|

|                      |                                               |
|----------------------|-----------------------------------------------|
| <b>Frame tubing:</b> | 3.2 cm (outside diameter), 1.5 mm (thickness) |
|----------------------|-----------------------------------------------|

|                             |                                  |
|-----------------------------|----------------------------------|
| <b>Swivel castor wheel:</b> | 2.40 × 7.2 cm (width × diameter) |
|-----------------------------|----------------------------------|

|                   |      |
|-------------------|------|
| <b>Tolerance:</b> | ± 4% |
|-------------------|------|

|                           |       |
|---------------------------|-------|
| <b>Carrying capacity:</b> | 40 kg |
|---------------------------|-------|

|                                 |     |
|---------------------------------|-----|
| <b>Knock-down construction:</b> | Yes |
|---------------------------------|-----|

OUR TRUSTED LEGACY IS YOUR SUCCESS



## HU-I-U0002156 SPECIFICATIONS

### SUPPLIED WITH:

- One (1) complete set of tools required for assembly.
- List of accessories and parts.
- Instructions for use, including cleaning and disinfection instructions, in English, French, and Spanish.

### PACKING, LABELLING, AND INSTRUCTIONS:

- One (1) unit per box.
- Unit weight: 11.90 kg (including packaging).
- Symbols in accordance with ISO 15223.