

OUR TRUSTED LEGACY IS YOUR SUCCESS



PRODUCT CATEGORY:
HOSPITAL FURNITURE

PRODUCT NAME:
HU-I-U0002157

PRODUCT DESCRIPTION:
INSTRUMENT TROLLEY



Area of Use: Operating Theatre

Intended Use: A mobile instrument table designed to support and position surgical instruments and supplies during procedures in the operating theatre.

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-I-U0002157 SPECIFICATIONS

TECHNICAL SPECIFICATION:

- Movable instrument table with two flat, non-removable shelves / trays.
- Heavy-duty carriage mounted on four swivel castors, of which two are fitted with brakes and two are fitted with anti-static wheels.
- Transfer bars connect all lower distal portions of the four base supports, providing maximum structural strength.
- Shelves have a smooth finish to allow easy cleaning.
- The top shelf has no raised edges, creating a completely flat work surface.
- The lower shelf has a fully flat surface and is fixed to the table frame directly above the swivel castors.

MATERIALS:

- Designed for high resistance to corrosion (tropical environments).
- Frame and trays: Austenitic stainless steel 18/10 (304 quality)
- Castor frame / bracket: Steel
- Castor brake: Total-lock type (wheel and swivel lock)
- Castor wheel: Single wheel, mould-on type, non-hooded (for easy maintenance), anti-static (for two wheels)
- Wheel bearings: Sealed bearings in both the swivel and the wheel
- Swivel mechanism: Ball-bearing
- All furniture components contain less than 100 parts per million (ppm) of the following chemical groups:
 - Urea formaldehyde.
 - Heavy metals, including mercury, cadmium, lead, and antimony.
 - Hexavalent chromium in plated finishes.
 - Stain- and non-stick treatments derived from Perfluorinated Compounds (PFCs), including Perfluorooctanoic Acid (PFOA).
 - Added antimicrobial treatments.

DIMENSIONS:

- Overall: 82 × 46 × 86.5 cm (L × W × H)
- Frame tubing: 3.2 cm (outside diameter), 1.5 mm (thickness)
- Shelf thickness: 1.5 mm
- Swivel castor wheel: 3.2 × 10 cm (width × diameter)
- Tolerance: ± 5%
- Carrying capacity: 100 kg
- Knock-down construction: Yes

SUPPLIED WITH:

- One (1) complete set of tools required for assembly.
- List of accessories and parts.
- Detailed, step-by-step assembly and safe-use instructions, text-free and pictorial (line drawings only).
- Instructions for use, including cleaning and disinfection instructions, in English, French, and Spanish.

OUR TRUSTED LEGACY IS YOUR SUCCESS

HU-I-U0002157 SPECIFICATIONS

PACKAGING, LABELLING, AND INSTRUCTIONS:
One (1) unit per box.
Unit weight: 23 kg (including packaging).
Unit volume: 0.06 m ³ (including packaging).
Symbols in accordance with ISO 15223.
Packaging made from recycled or recyclable materials.