

OUR TRUSTED LEGACY IS YOUR SUCCESS



PRODUCT CATEGORY:
HOSPITAL FURNITURE

PRODUCT NAME:
HU-I-U0002159

PRODUCT DESCRIPTION:
HOSPITAL BED ICU



Area of Use: All Areas

Intended Use: A manually adjustable hospital bed intended to provide support and positioning for patients during rest, sleep, treatment, and recovery in healthcare settings.

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-I-U0002159 SPECIFICATIONS

TECHNICAL SPECIFICATION:

Hospital bed, intensive care, four sections.

Mounted on four heavy-duty swivel castors, two of which are fitted with brakes.

Transfer bars connect the lower distal portions of the two foot-end legs and the two head-end legs, providing maximum structural strength.

Protective bumpers are fitted at all four corners.

Bed ends are finished with panels.

Four-section mattress platform made of epoxy-painted steel mesh with side supports to immobilize the mattress.

Mattress cover is removable via a side zipper.

Manually operated crank allows adjustment of the backrest (up to 80°), leg section, and foot section (up to 35°).

The crank handle folds away underneath the bed when not in use.

Adjustable and removable folding side rails are provided.

MATERIALS:

Designed for high resistance to corrosion (tropical environments).

Frame:	Epoxy-coated tubular steel
---------------	----------------------------

Mattress:	High-density polyurethane foam, density 27–30 kg/m ³
------------------	---

Cover:	Flexible plastic, highly tear-resistant, anti-static, flame-retardant, non-absorbing, waterproof, and cleanable with hospital-grade disinfectants
---------------	---

Castor frame / bracket:	Steel
--------------------------------	-------

Castor wheel:	Single wheel, mould-on type, non-hooded (for easy maintenance)
----------------------	--

Wheel bearings:	Sealed bearings in both the swivel and the wheel
------------------------	--

Swivel mechanism:	Ball-bearing
--------------------------	--------------

All furniture components contain less than 100 parts per million (ppm) of the following chemical groups:

Urea formaldehyde.

Heavy metals, including mercury, cadmium, lead, and antimony.

Hexavalent chromium in plated finishes.

Stain- and non-stick treatments derived from Perfluorinated Compounds (PFCs), including Perfluorooctanoic Acid (PFOA).

Added antimicrobial treatments.

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-I-U0002159 SPECIFICATIONS

DIMENSIONS:

Sleeping surface:	200 × 85 cm (L × W)
Height of surface (without mattress):	50 cm (fixed)
Mattress thickness:	12 cm
Bed frame tubing:	5 × 3 cm (H × W), 2 mm (thickness)
Leg frame tubing:	3.2 cm diameter, 2 mm (thickness)
Swivel castor wheels:	3.2 × 12.4 cm (width × diameter)
Tolerance:	±4%
Castors with brakes: 2	2
Carrying capacity:	165 kg
Knock-down construction:	Yes

SUPPLIED WITH:

One (1) complete set of tools required for assembly.
One (1) matching mattress with cover.
Two (2) removable folding side rails.
List of accessories and parts.
Detailed, step-by-step assembly and safe-use instructions, text-free and pictorial (line drawings only).
Instructions for use, including cleaning and disinfection instructions, in English, French, and Spanish.

PACKAGING, LABELLING, AND INSTRUCTIONS:

One (1) unit per box.
Unit weight: 77.1 kg (including packaging).
Unit volume: 0.51 m ³ (including packaging).
Symbols in accordance with ISO 15223.
Packaging made from recycled or recyclable materials.