

OUR TRUSTED LEGACY IS YOUR SUCCESS



PRODUCT CATEGORY:
HOSPITAL FURNITURE

PRODUCT NAME:
HU-I-U0002161

PRODUCT DESCRIPTION:
KICK BUCKET



Area of Use: All Areas

Intended Use: A kick bucket is a circular frame mounted on wheels, with a bucket which is used for the receiving of medical waste. It could easily be moved around by a persons foot, for example in an operating room.

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-I-U0002161 SPECIFICATIONS

TECHNICAL SPECIFICATION:

Waste bucket with handle, mounted on a heavy-duty mobile holder frame.

The bucket handle pivots and is fully integrated through robust hinges.

Wide, heavy base mounted on five swivel castors, all fitted with anti-static wheels.

The bucket can be smoothly placed into and removed from the holder frame.

The holder frame is fitted with a 360-degree rubber bumper for protection.

MATERIALS:

Designed for high resistance to corrosion (tropical environments).

Frame and bucket:	Austenitic stainless steel 18/10 (304 quality)
-------------------	--

Bumper:	Rubber
---------	--------

Castor frame / bracket:	Steel
-------------------------	-------

Castor wheel:	Single wheel, mould-on type, non-hooded (for easy maintenance), anti-static (all five)
---------------	--

Swivel mechanism:	Ball-bearing
-------------------	--------------

All furniture components contain less than 100 parts per million (ppm) of the following chemical groups, as applicable:

Urea formaldehyde.

Heavy metals, including mercury, cadmium, lead, and antimony.

Hexavalent chromium in plated finishes.

Stain- and non-stick treatments derived from Perfluorinated Compounds (PFCs), including Perfluorooctanoic Acid (PFOA).

Added antimicrobial treatments.

DIMENSIONS:

Overall:	38 × 38 cm (diameter × height)
----------	--------------------------------

Swivel castor wheel:	2 × 5 cm (width × diameter)
----------------------	-----------------------------

Tolerance:	±5%
------------	-----

Bumper height:	3 cm
----------------	------

Bucket capacity:	11.7 litres
------------------	-------------

Tolerance:	±5%
------------	-----

Knock-down construction:	No
--------------------------	----

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-I-U0002161 SPECIFICATIONS

SUPPLIED WITH:

- One (1) complete set of tools required for assembly.
- List of accessories and parts.
- Detailed, step-by-step assembly and safe-use instructions, text-free and pictorial (line drawings only).
- Instructions for use, including cleaning and disinfection instructions, in English, French, and Spanish.

PACKAGING, LABELLING, AND INSTRUCTIONS:

- One (1) unit per box.
- Unit weight: 7.1 kg (including packaging).
- Unit volume: 0.085 m³ (including packaging).
- Symbols in accordance with ISO 15223.
- Packaging made from recycled or recyclable materials.