

OUR TRUSTED LEGACY IS YOUR SUCCESS



PRODUCT CATEGORY:
HOSPITAL FURNITURE

PRODUCT NAME:
HU-I-U00021562

PRODUCT DESCRIPTION:
EMERGENCY TROLLEY



Area of Use: All Areas

Intended Use: The mobile trolley/cart is intended for the storage and transport of medical devices, medications, and related supplies within healthcare facilities. It is equipped with drawers for secure storage and a top surface to support and organise items required for clinical procedures. The trolley is mobile, allowing convenient positioning at the point of care to assist healthcare professionals during medical interventions.

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-I-U0002162 SPECIFICATIONS

TECHNICAL SPECIFICATION:

Emergency response trolley with work surface and storage space.

Heavy-duty carriage mounted on four swivel castors, two fitted with brakes and all four fitted with anti-static wheels.

The work surface has raised edges along one length and both widths and is finished with an anti-slip layer.

One long side is fitted with four drawers arranged side-by-side for storage of medicines, consumables, and small equipment.

The front of each drawer is fitted with a pre-fixed holder for content identification strips.

Drawers have recessed finger grips, leaving the front entirely flat.

One central locking system secures all drawers.

Each trolley is supplied with a unique key. The key cannot be broken while inside the lock (fold-away or similar design).

The inside of the drawers is customizable with organiser dividers.

Integrated fittings with lids are provided for a waste container (total capacity: 8 litres) and a sharps container (total capacity: 2.5 litres).

A flat worktop extension slides in and out from underneath the main work surface, extending the workspace.

One short side is fitted with a push-bar handle.

Both braked castors are mounted on the same side as the push-bar handle.

Protective bumpers are fitted at all four corners.

MATERIALS:

Designed for high resistance to corrosion (tropical environments).

Frame, side panels, base, and drawers: Epoxy-coated steel plate

Push handle: Austenitic stainless steel 18/10

Worktop and extension: ABS

Castor frame / bracket: Steel

Castor brake: Total-lock type (wheel and swivel lock)

Castor wheel: Single wheel, mould-on type, non-hooded (for easy maintenance), anti-static (all four)

Wheel bearings: Sealed bearings in both the swivel and the wheel

Swivel mechanism: Ball-bearing

All furniture components contain less than 100 parts per million (ppm) of the following chemical groups:

Urea formaldehyde.

Heavy metals, including mercury, cadmium, lead, and antimony.

Hexavalent chromium in plated finishes.

Stain and non-stick treatments derived from Perfluorinated Compounds (PFCs), including Perfluorooctanoic Acid (PFOA).

Added antimicrobial treatments.

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-I-U0002162 SPECIFICATIONS

DIMENSIONS:	
Trolley, overall:	81 × 60 × 100 cm (L × W × H)
Frame, drawers, and panels thickness:	1 mm
Raised edges, work surface:	4.8 cm (height)
Worktop extension:	40 × 36 cm (L × W)
Upper drawers height:	10 cm
Second drawers height:	15 cm
Third drawers height:	15 cm
Base drawers height:	30 cm
Swivel castor wheel:	3.2 × 12.5 cm (width × diameter)
Tolerance:	±4%
Carrying capacity:	120 kg
Knock-down construction:	No

SUPPLIED WITH:
One (1) complete set of tools required for assembly.
One (1) set of organisers for each drawer.
One (1) set of two keys, unique per trolley.
List of accessories and parts.
Detailed, step-by-step assembly and safe-use instructions, text-free and pictorial (line drawings only).

PACKAGING, LABELLING, AND INSTRUCTIONS:
One (1) unit per box.
Unit weight: 65 kg (including packaging).
Unit volume: 0.59 m ³ (including packaging).
Symbols in accordance with ISO 15223.
Packaging made from recycled or recyclable materials.