

OUR TRUSTED LEGACY IS YOUR SUCCESS



PRODUCT CATEGORY:
HOSPITAL FURNITURE

PRODUCT NAME:
HU-I-U0002165

PRODUCT DESCRIPTION:
BEDSIDE CABINET



Area of Use: All Areas

Intended Use: A bedside cabinet provides flexible and convenient storage for the patients' belongs, as well as being a fixed bedside tabletop.

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-I-U0002165 SPECIFICATIONS

TECHNICAL SPECIFICATION:

Movable basic patient bedside cabinet.
Mounted on 4 swivel castors, two with brake.
Lower part: storage compartment with one fixed shelf, closed with door.
Upper part: drawer, lockable with key.
Key is unique for each single cabinet, key-bow folds away from key-blade avoiding breakage when inside lock.
One side rail handle allows for easy repositioning.
Both castors with brake are mounted at the push-bar handle side.
Drawers and cabinet have recessed finger grip-hold, leaving front entirely flat.

MATERIALS:

High resistance to corrosion (tropical environment).	
Frame, side panels, base, top, door and shelf:	Epoxy coated plate steel.
Top:	Has smooth finishing allowing for easy cleaning.
Caster frame/bracket:	Nylon.
Caster brake:	Wheel lock.
Caster wheel:	Double wheels, mold-on type, non-hooded (for easy maintenance).
All furniture components contain less than 100 parts per million of the below chemical groups:	
Urea formaldehyde.	
Heavy metals including mercury, cadmium, lead, antimony.	
Hexavalent chromium in plated finishes.	
Stain and non-stick treatments derived from Perfluorinated Compounds (PFCs), including Perfluorooctanoic Acid (PFOA).	
Added antimicrobial treatments.	

DIMENSIONS:

Cabinet, overall:	40 x 40 x 80 cm (l x w x h)
Frame, drawer, door and panels:	0,8 mm thickness
Door, opening angle:	120 degree
Swivel castor wheel:	(4) x (5.2) cm (w x diameter) (Plus and minus %4)
Carrying capacity:	30 kg
Knockdown construction:	No

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-I-U0002165 SPECIFICATIONS

SUPPLIED WITH:

- 1 x complete set of tools required for assembly.
- 1 x set of 2 keys, unique per cabinet .
- 1 x set of organisers for each drawer .
- List of accessories and parts .
- Detailed step-by-step instructions for assembly and safe use, text-free pictorial based (i.e. line-drawings only).
- Instructions for use including cleaning instructions in English, French and Spanish.

PACKING, LABELLING, INSTRUCTIONS:

- One (1) unit per box.
- Unit Weight in 15.6 Kg (including its packaging).
- Unit Volume in 0.16 m3 (including its packaging).
- Symbols in accordance with ISO 15223.
- The packaging is made of recycled or recyclable materials.