

OUR TRUSTED LEGACY IS YOUR SUCCESS



PRODUCT CATEGORY:
HOSPITAL FURNITURE

PRODUCT NAME:
HU-I-U0002166

PRODUCT DESCRIPTION:
DOUBLE DOOR INSTRUMENT STORAGE CABINET



Area of Use: Normal Ward to ICU.

Intended Use: The cabinet is intended for the secure storage of medical instruments and related supplies within healthcare facilities. It is designed to protect contents from contamination, damage, and unauthorised access, while ensuring organised and efficient storage.

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-I-U0002166 SPECIFICATIONS

TECHNICAL SPECIFICATION:

Instrument storage cabinet, double door.

Mounted on four sturdy supports, all fitted with height-adjustable feet.

Recessed safety glass in the doors allows viewing of the cabinet contents.

The safety glass is supported by the door frame on both sides, leaving no exposed glass edges.

Reinforced internal fixtures allow height adjustment of the four shelves.

Doors are triple-hinged and fitted with handles and locks.

Two-point locking system secures the doors at both the top and bottom of the frame.

Each cabinet is supplied with a unique key. The key bow folds away from the key blade to prevent breakage while inside the lock.

Clearance underneath the cabinet allows for easy cleaning.

MATERIALS:

Designed for high resistance to corrosion (tropical environments).

Adjustable feet: Nylon

Frame, side panels, base, top, and shelves: Epoxy-coated steel plate

Doors: Hardened safety glass with ground edges, fitted with key lock

All furniture components contain less than 100 parts per million (ppm) of the following chemical groups:

Urea formaldehyde.

Heavy metals, including mercury, cadmium, lead, and antimony.

Hexavalent chromium in plated finishes.

Stain- and non-stick treatments derived from Perfluorinated Compounds (PFCs), including Perfluorooctanoic Acid (PFOA).

Added antimicrobial treatments.

DIMENSIONS:

Cabinet, overall: 80 × 40 × 190 cm (L × W × H)

Clearance underneath: 10 cm

Safety glass thickness: 3 mm

Front glass area: 70%

Frame, door, and panels thickness: 1 mm

Door opening angle: 180°

Carrying capacity per shelf: 25 kg

Knock-down construction: Yes

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-I-U0002166 SPECIFICATIONS

SUPPLIED WITH:
One (1) complete set of tools required for assembly.
One (1) set of four shelves.
One (1) set of two keys, unique per cabinet.
List of accessories and parts.
Detailed, step-by-step assembly and safe-use instructions, text-free and pictorial (line drawings only).
Instructions for use, including cleaning instructions, in English, French, and Spanish.
PACKAGING, LABELLING, AND INSTRUCTIONS:
One (1) unit per box.
Unit weight: 68 kg (including packaging).
Unit volume: 0.27 m ³ (including packaging).
Symbols in accordance with ISO 15223.
· Packaging made from recycled or recyclable materials