

OUR TRUSTED LEGACY IS YOUR SUCCESS



PRODUCT CATEGORY:
HOSPITAL FURNITURE

PRODUCT NAME:
HU-I-U0002168

PRODUCT DESCRIPTION:
SINGLE BOWL STAND



Area of Use: Operating Theatres, Procedure Rooms

Intended Use: The mobile stand with single bowl is intended to support and hold a bowl for use during surgical and clinical procedures. It is designed to contain fluids, medical waste, or surgical materials as required, and to provide convenient access for healthcare professionals. The stand is mobile, allowing easy positioning close to the patient within operating theatres and procedure rooms.

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-I-U0002168 SPECIFICATIONS

TECHNICAL SPECIFICATION:

Movable, fixed-height stand designed to hold a single bowl.
Mounted on a heavy-duty carriage base fitted with five swivel castors for smooth mobility and stability.
Lower support column is deeply and securely fixed into the carriage base to ensure structural strength.
Upper section of the support column is fitted with a ring structure designed to securely hold a compatible bowl.
Bowl capacity: 3 litres.

MATERIALS:

Support column, ring structure, and bowl:	Austenitic stainless steel 18/10
Carriage base:	Single moulded polypropylene unit
Castor frame/bracket:	Polypropylene
Castor wheel:	Twin-wheel, solid type, non-hooded (for easy maintenance), polypropylene
Wheel bearing:	Polypropylene
Castor fixation into carriage base:	Stainless steel

All furniture components contain less than 100 parts per million (ppm) of the following chemical groups:

Urea formaldehyde

Heavy metals, including mercury, cadmium, lead, and antimony

Hexavalent chromium in plated finishes

Stain and non-stick treatments derived from Perfluorinated Compounds (PFCs), including Perfluorooctanoic Acid (PFOA)

Added antimicrobial treatments

DIMENSIONS:

Stand height:	100 cm
Bowl holder diameter:	30 cm
Carriage base diameter:	60 cm
Support column:	3.2 cm outer diameter, 1.5 mm thickness
Swivel castor twin-wheel:	4 × 5 cm (W × diameter) (±3%)
Carrying capacity:	10 kg
Knock-down construction:	Yes

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-I-U0002168 SPECIFICATIONS

SUPPLIED WITH:

- One (1) complete set of tools required for assembly.
- One (1) stainless steel bowl, 3 L capacity.
- List of accessories and parts.
- Detailed step-by-step assembly and safe-use instructions, text-free pictorial format (line drawings only).
- Instructions for use, including cleaning and disinfection instructions, in English, French, and Spanish.

PACKAGING, LABELLING, AND INSTRUCTIONS:

- One (1) unit per box.
- Unit weight: 5.6 kg (including packaging).
- Unit volume: 0.07 m³ (including packaging).
- Symbols in accordance with ISO 15223.
- Packaging made from recycled or recyclable materials.