

OUR TRUSTED LEGACY IS YOUR SUCCESS



PRODUCT CATEGORY:
HOSPITAL FURNITURE

PRODUCT NAME:
HU-I-U0002169

PRODUCT DESCRIPTION:
FOOTSTOOL, TWO STEPS.



Area of Use: All Areas

Intended Use: The two-step footstool is intended to provide a stable and elevated platform to assist hospital staff and patients in safely reaching higher surfaces or accessing examination tables, beds, or equipment. It supports safe mounting and dismounting and is suitable for use in clinical and healthcare environments.

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-I-U0002169 SPECIFICATIONS

TECHNICAL SPECIFICATION:

- Sturdy two-step footstool.
- Mounted on robust supporting legs, widely spaced to provide maximum stability.
- Steps feature an embossed anti-slip pattern, and the feet are fitted with heavy-duty anti-slip nylon ends.
- Transfer bars connect all lower distal portions of the supports, providing maximum structural strength.

MATERIALS:

- High resistance to corrosion (tropical environment).
- Anti-slip feet: Rubber.
- Frame: Epoxy-coated tubular steel.
- Steps: Heavy-duty solid plastic.
- All furniture components contain less than 100 parts per million of the following chemical groups:
 - Urea formaldehyde.
 - Heavy metals including mercury, cadmium, lead, and antimony.
 - Hexavalent chromium in plated finishes.
 - Stain and non-stick treatments derived from perfluorinated compounds (PFCs), including perfluorooctanoic acid (PFOA).
 - Added antimicrobial treatments.

DIMENSIONS:

Footstool, overall:	62 x 38 x 42 cm (L x W x H)
Frame:	2.5 cm (outside, across), 1.5 mm (thickness)
Carrying capacity:	100 kg
Knock-down construction:	Yes

SUPPLIED WITH:

- 1 x Complete set of tools required for assembly.
- List of accessories and parts.
- Detailed step-by-step instructions for assembly and safe use, text-free pictorial-based (i.e. line drawings only).
- Instructions for use, including cleaning instructions, in English, French, and Spanish.

PACKING, LABELLING, INSTRUCTIONS:

- One (1) unit per box
- Unit weight: 8 kg (including packaging)
- Unit volume: 0.032 m³ (including packaging)
- Symbols in accordance with ISO 15223
- Packaging made of recycled or recyclable materials