

OUR TRUSTED LEGACY IS YOUR SUCCESS



PRODUCT CATEGORY:
HOSPITAL FURNITURE

PRODUCT NAME:
HU-I-U0002170

PRODUCT DESCRIPTION:
PATIENT STRETCHER WITH RAILS



Area of Use: All Areas

Intended Use: Manually operated patient stretcher designed for the safe transport, transfer, and temporary positioning of patients within healthcare facilities, including hospitals, emergency departments, and clinical areas. Equipped with side rails to enhance patient safety during movement and transfer.

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-I-U0002170 SPECIFICATIONS

TECHNICAL SPECIFICATION:	
Stretcher designed for safe transport of patients within healthcare facilities, two-section type.	
Mounted on a heavy-duty carriage fitted with four swivel castors, two equipped with brakes. All castors are fitted with anti-static wheels.	
Both sections are fitted with fixed, padded upholstery.	
Backrest section is adjustable via a secure pawl and gear ratchet mechanism, ensuring safe operation for both patient and caregiver.	
When fully extended, both sections align to form a completely flat patient surface.	
Transfer bars connect the lower portions of all four castors to provide maximum structural strength and stability.	
Base of the stretcher is fitted with a large, securely fixed mesh utility shelf for storage of medical items during transport. Shelf remains stable without shaking or vibration during movement.	
Equipped with fold-away side rails for patient safety during transport.	
Fitted with protective bumpers at all four corners to prevent damage to walls and equipment.	
Head-end fitted with removable, height-adjustable IV pole. Height adjustment is secured using a robust clamp with heavy-duty adjustment knob.	
Fixings for side rails and IV pole are solid steel and permanently welded to the stretcher frame.	
MATERIALS:	
Designed for high resistance to corrosion in tropical environments.	
Frame:	Epoxy-coated tubular steel
Upholstery padding:	High-density polyurethane foam, density 27–30 kg/m ³
Upholstery cover:	Medical-grade plastic material, flexible, highly tear-resistant, anti-static, flame-retardant, non-absorbing, waterproof, and compatible with hospital-grade cleaning and disinfection products.
Castor frame/bracket:	Steel
Castor brake:	Total-lock type (simultaneous wheel and swivel lock)
Castor wheel:	Single-wheel, mould-on type, non-hooded, anti-static
Wheel bearing:	Sealed bearings in both swivel and wheel
Swivel mechanism:	Ball-bearing type
All furniture components contain less than 100 parts per million (ppm) of the following chemical groups:	
Urea formaldehyde.	
Heavy metals including mercury, cadmium, lead, and antimony.	
Hexavalent chromium in plated finishes.	
Perfluorinated Compounds (PFCs), including Perfluorooctanoic Acid (PFOA).	
Added antimicrobial treatments.	

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-I-U0002170 SPECIFICATIONS

DIMENSIONS:	
Overall dimensions (including upholstery):	180 × 55 × 80 cm (L × W × H)
Fold-away side rails:	80 × 35 cm (L × H)
Frame tubing:	3.2 cm outer diameter, 2 mm thickness
Castor wheels:	3.2 × 12.5 cm (W × diameter) (±5%)
Upholstery thickness:	5 cm
Maximum safe working load:	160 kg
Knock-down construction:	Yes
SUPPLIED WITH:	
One (1) complete set of tools required for assembly	
One (1) mesh utility shelf	
Two (2) fold-away side rails	
One (1) IV pole	
List of accessories and parts	
Step-by-step assembly and safe-use instructions in text-free pictorial format (line drawings only)	
Instructions for use, including cleaning and disinfection instructions, in English, French, and Spanish	
PACKAGING, LABELLING, AND INSTRUCTIONS:	
One (1) unit per box	
Unit weight: 47 kg (including packaging)	
Unit volume: 0.256 m³ (including packaging)	
Symbols in accordance with ISO 15223	
Packaging made from recycled or recyclable materials	