

OUR TRUSTED LEGACY IS YOUR SUCCESS



PRODUCT CATEGORY:
HOSPITAL FURNITURE

PRODUCT NAME:
HU-I-U0002171

PRODUCT DESCRIPTION:
EXAMINATION TABLE



Area of Use: All Areas

Intended Use: An examination table is used to support patients during medical examinations. During these exams, doctors, clinics and hospitals use an adjusting mechanism to manipulate and position the table to allow patient support, closer examination of the patient, and to allow the patient to move on and off the table safely.

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-I-U0002171 SPECIFICATIONS

TECHNICAL SPECIFICATION:

- Examination table, two-section type, designed for patient examination and minor medical procedures.
- Mounted on four sturdy supports, each fitted with height-adjustable feet to ensure stability on uneven surfaces.
- Both sections are fitted with fixed, padded upholstery.
- Backrest section is adjustable via a secure pawl and gear ratchet mechanism, ensuring safe and reliable positioning for both patient and caregiver.
- When fully extended, both sections align to form a completely flat examination surface.
- Transfer bars connect the lower portions of all four supports, providing maximum structural strength and stability.

MATERIALS:

- Designed for high resistance to corrosion in tropical environments.
- Frame:** Epoxy-coated tubular steel.
- Adjustable feet:** Nylon
- Upholstery padding:** High-density polyurethane foam, density 27–30 kg/m³.
- Upholstery cover:** Medical-grade plastic material, flexible, highly tear-resistant, anti-static, flame-retardant, non-absorbing, waterproof, and compatible with hospital-grade cleaning and disinfection products.
- All furniture components contain less than 100 parts per million (ppm) of the following chemical groups:
 - Urea formaldehyde.
 - Heavy metals including mercury, cadmium, lead, and antimony.
 - Hexavalent chromium in plated finishes.
 - Perfluorinated Compounds (PFCs), including Perfluorooctanoic Acid (PFOA).
 - Added antimicrobial treatments.

DIMENSIONS:

- | | |
|--|---|
| Overall dimensions (including upholstery): | 190 × 55 × 80 cm (L × W × H) |
| Frame tubing: | 3.2 cm outer diameter, minimum 1.5 mm thickness |
| Upholstery thickness: | 5 cm |
| Maximum safe working load: | 160 kg |
| Knock-down construction: | Yes |

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-I-U0002171 SPECIFICATIONS

SUPPLIED WITH:

- One (1) complete set of tools required for assembly.
- List of accessories and parts.
- Step-by-step assembly and safe-use instructions in text-free pictorial format (line drawings only).
- Instructions for use, including cleaning and disinfection instructions, in English, French, and Spanish.

PACKAGING, LABELLING, AND INSTRUCTIONS:

- One (1) unit per box
- Unit weight: 30.2 kg (including packaging).
- Unit volume: 0.193 m³ (including packaging).
- Symbols in accordance with ISO 15223.
- Packaging made from recycled or recyclable materials.