

OUR TRUSTED LEGACY IS YOUR SUCCESS



PRODUCT CATEGORY:
HOSPITAL FURNITURE

PRODUCT NAME:
HU-I-U0002174

PRODUCT DESCRIPTION:
DELIVERY TABLE



Area of Use: Maternity Wards / Delivery Rooms

Intended Use: A manually operated gynecology and delivery table designed to support and position a woman safely and comfortably during labour, delivery, and gynaecological examinations and procedures related to pregnancy and childbirth in healthcare facilities.

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-I-U0002174 SPECIFICATIONS

TECHNICAL SPECIFICATIONS:

- Gynecological examination and delivery table, three sections.
- Mounted on four sturdy supports, all fitted with height-adjustable feet.
- All sections are fitted with non-removable padded upholstery.
- Robust mechanics allow for manual repositioning of all sections between gynecological and obstetric use.
- Transfer bars connect all lower distal portions of the four supports, providing maximum structural strength.

Back Section:

- Adjustable via secured pawl and gear ratchet, ensuring safety for both patient and operator.

Pelvic Section:

- Can be tilted to the Trendelenburg position using a pawl and gear ratchet.
- Sides of this section are fitted with vertical handgrips, length 20 cm.
- Soft-touch plastic knee crutches are height and width adjustable, secured with robust clamps and heavy knobs.
- Fixing of the crutch holders is solid steel and welded to the frame of the bed.
- Pelvic section integrates a support for a slide-out basin tray.

Leg Section:

- Recesses entirely downwards to 90 degrees.
- When elevated and fully extended, all sections align to form a perfectly flat surface.

MATERIALS:

- High resistance to corrosion (tropical environment).
- Frame:** Epoxy-coated tubular steel.
- Adjustable feet:** Nylon.
- Sliders / fixtures for the knee crutches:** Tubular steel, welded to the bed frame.
- Padded upholstery:** High-density polyurethane foam, density 27–30 kg/m³.
- Cover:** Plastic, flexible, highly tear-resistant, anti-static, flame-retardant, non-absorbing, waterproof, and cleanable with hospital-grade disinfection products.
- Basin tray:** Austenitic stainless steel, 18/10.
- All furniture components contain less than 100 parts per million of the following chemical groups:
 - Urea formaldehyde.
 - Heavy metals including mercury, cadmium, lead, and antimony.
 - Hexavalent chromium in plated finishes.
 - Stain and non-stick treatments derived from perfluorinated compounds (PFCs), including perfluorooctanoic acid (PFOA).
 - Added antimicrobial treatments.

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-I-U0002174 SPECIFICATIONS

DIMENSIONS:	
All three sections extended, including padding:	180 x 70 x 80 cm (L x W x H)
Frame:	3.2 cm (outside, across), 2 mm (thickness)
Padded upholstery:	5 cm (H)
Carrying capacity:	150 kg
Slide-out basin tray capacity:	4 L
Knock-down construction:	Yes

PACKAGING, LABELLING, INSTRUCTIONS:	
One (1) unit per box.	
Unit weight: 46 kg (including packaging).	
Unit volume: 0.23 m ³ (including packaging).	
Symbols in accordance with ISO 15223.	
Packaging made of recycled or recyclable materials.	

SUPPLIED WITH:	
1 x Complete set of tools required for assembly.	
2 x Leg holders with canvas straps, adjustable height and width.	
2 x Knee crutches, adjustable height and width.	
1 x Basin tray.	
List of accessories and parts.	
Detailed step-by-step instructions for assembly and safe use, text-free pictorial-based (i.e. line drawings only).	
Instructions for use, including cleaning and disinfection instructions, in English, French, and Spanish.	