

OUR TRUSTED LEGACY IS YOUR SUCCESS



PRODUCT CATEGORY:
HOSPITAL FURNITURE

PRODUCT NAME:
HU-I-U0002169

PRODUCT DESCRIPTION:
BLOOD DONATION CHAIR



Area of Use: Procedure Rooms

Intended Use: A manually operated reclining chair intended to position and support blood donors during blood collection procedures, ensuring patient comfort, stability, and appropriate clinical access for healthcare personnel.

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-I-U0002169 SPECIFICATIONS

TECHNICAL SPECIFICATION:

Blood donation chair.

Mounted on four sturdy supports, all fitted with height-adjustable feet.

Fixed height.

All sections are fitted with non-removable padded upholstery.

Robust mechanisms allow for manual repositioning of all sections.

Self-locking adjustable armrest.

MATERIALS:

High resistance to corrosion (tropical environment).

Frame: Epoxy-coated steel.

Adjustable feet: Nylon.

Padded upholstery: High-density polyurethane foam, density 27–30 kg/m³.

Cover: Plastic, flexible, highly tear-resistant, anti-static, flame-retardant, non-absorbing, waterproof, and cleanable with hospital-grade disinfection products.

All furniture components contain less than 100 parts per million of the following chemical groups:

Urea formaldehyde.

Heavy metals including mercury, cadmium, lead, and antimony.

Hexavalent chromium in plated finishes.

Stain and non-stick treatments derived from Perfluorinated Compounds (PFCs), including Perfluorooctanoic Acid (PFOA).

Added antimicrobial treatments.

DIMENSIONS:

Width: 55 cm

Height: 65 cm

Padded upholstery thickness: 5.5 cm

Length (fully extended position): 187 cm

Carrying capacity: 150 kg

Knockdown construction: Yes

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-I-U0002169 SPECIFICATIONS

SUPPLIED WITH:
1 × Complete set of tools required for assembly.
List of accessories and parts.
Detailed step-by-step instructions for assembly and safe use, text-free pictorial-based (i.e. line drawings only).
Instructions for use, including cleaning and disinfection instructions, in English, French, and Spanish.
PACKAGING, LABELLING, AND INSTRUCTIONS:
Unit weight: 30.5 kg
Unit volume: 0.16 m³
Symbols in accordance with ISO 15223
Packaging is made of recycled or recyclable materials.