


OUR TRUSTED LEGACY IS YOUR SUCCESS



PRODUCT CATEGORY:
HOSPITAL FURNITURE

PRODUCT NAME:
HU-I-U0101015

PRODUCT DESCRIPTION:
PATIENT PRIVACY SCREEN



Area of Use: All Areas

Intended Use: The hospital bedscreen is intended to provide visual privacy for patients during medical examinations, treatments, or care procedures. It may also be used as a room divider to separate patient areas within healthcare facilities. The screen is mobile, allowing easy positioning around patient beds or examination tables as required.

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-I-U0101015 SPECIFICATIONS

TECHNICAL SPECIFICATIONS:

Hospital bedscreen, three sections with curtains.

Bedscreen mounted on swivel castors and positioned in a broad stance for good stability.

Curtains are integrated centrally into the frame, leaving no gaps near the hinges.

MATERIALS:

High resistance to corrosion (tropical environment).

Frame: Epoxy-coated plated tubular steel.

Curtains: Plastic, non-transparent, tear-resistant, flame-retardant, disinfectant- and liquid-proof, washable, light-coloured.

Caster frame / bracket: Nylon.

All furniture components contain less than 100 parts per million of the following chemical groups:

Urea formaldehyde.

Heavy metals including mercury, cadmium, lead, and antimony.

Hexavalent chromium in plated finishes.

Stain and non-stick treatments derived from perfluorinated compounds (PFCs), including perfluorooctanoic acid (PFOA).

Added antimicrobial treatments.

DIMENSIONS:

Bedscreen: Three panels, each 70 x 170 cm (W x H)

Swivel castors: 5 cm ($\pm 4\%$)

Knock-down construction: Yes

SUPPLIED WITH:

1 x Complete set of tools required for assembly

1 x Set of plastic curtains

List of accessories and parts

Detailed step-by-step instructions for assembly and safe use, text-free pictorial-based (i.e. line drawings only)

Instructions for use, including cleaning instructions, in English, French, and Spanish

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-I-U0101015 SPECIFICATIONS

PACKAGING, LABELLING, INSTRUCTIONS:
One (1) unit per box.
Unit weight: 27.3 kg (including packaging).
Unit volume: 0.195 m³ (including packaging).
Meets requirements for palletisation, packaging, and labelling as specified here:
Symbols in accordance with ISO 15223.
Packaging made of recycled or recyclable materials.