


OUR TRUSTED LEGACY IS YOUR SUCCESS



PRODUCT CATEGORY:
HOSPITAL FURNITURE

PRODUCT NAME:
HU-I-U0120021

PRODUCT DESCRIPTION:
BASIC MEDICAL COT WITH IV POLE



Area of Use: All Areas

Intended Use: The manual foldable medical bed is intended to provide a safe and stable surface for patient rest, sleep, and recovery in hospital and healthcare environments. The bed is manually operated and may be folded for storage and transport when not in use.

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-I-U0120021 SPECIFICATIONS

TECHNICAL SPECIFICATIONS:

Basic medical cot/bed.

Two-section platform.

Backrest elevation with at least three positions (20, 30, and 45 degrees).

Frame is fitted with a sleeve to receive the IV pole.

IV pole, adjustable height range of 120 to 200 cm, with screw height adjustment.

IV pole fitted with double hooks.

MATERIALS:

High resistance to corrosion (tropical environment).

Frame and mesh:	Epoxy-coated steel.
------------------------	---------------------

Mattress:	Thickness 6.5 cm, medium-density polyurethane foam, density 26–30 kg/m ³ .
------------------	---

Cover:	Plastic, flexible, highly tear-resistant, anti-static, flame-retardant, non-absorbing, waterproof, and cleanable with hospital-grade disinfection products.
---------------	---

All furniture components contain less than 100 parts per million of the following chemical groups:

Urea formaldehyde.

Heavy metals including mercury, cadmium, lead, and antimony.

Hexavalent chromium in plated finishes.

Stain and non-stick treatments derived from perfluorinated compounds (PFCs), including perfluorooctanoic acid (PFOA).

Added antimicrobial treatments.

DIMENSIONS:

Open:	185 x 64 x 50 cm (L x W x H)
--------------	------------------------------

Folded:	90 x 64 x 30 cm (L x W x H)
----------------	-----------------------------

Carrying capacity:	120 kg
---------------------------	--------

Knock-down construction:	Yes
---------------------------------	-----

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-I-U0120021 SPECIFICATIONS

SUPPLIED WITH:

1 x Complete set of tools required for assembly
List of accessories and parts
Detailed step-by-step instructions for assembly and safe use, text-free pictorial-based (i.e. line drawings only)
1 x Fixed-height IV pole (multiple-sectioned for compact packaging)
1 x Mattress with cover
Instructions for use, including cleaning and disinfection instructions, in English, French, and Spanish

PACKAGING, LABELLING, INSTRUCTIONS:

One (1) unit per box.
Unit weight: 32.5 kg (including packaging).
Unit volume: 0.21 m ³ (including packaging).
Symbols in accordance with ISO 15223.
Packaging made of recycled or recyclable materials.