

OUR TRUSTED LEGACY IS YOUR SUCCESS



PRODUCT CATEGORY:  
**MEDICAL ACCESSORIES**

PRODUCT NAME:  
**HU-KBD**

PRODUCT DESCRIPTION:  
**WIRELESS KEYBOARD (MEDICAL GRADE)**



**Area of Use:** ALL AREAS

**Intended Use:** The wireless medical-grade keyboard is designed for use in clinical and healthcare environments to facilitate the safe and efficient input of data into medical devices, computers, or information systems. Its sealed, easy-to-clean design helps maintain infection control protocols by allowing thorough disinfection between uses, reducing the risk of cross-contamination. This device is intended for use by trained healthcare professionals in areas such as hospitals, clinics, laboratories, and other medical settings.

OUR TRUSTED LEGACY IS YOUR SUCCESS



## HU-KBD FEATURES:

- Wireless Silicone Keyboard & Touchpad
- IP67 Waterproof & Dustproof Protection Level
- 104 Keyboard & Function Keys
- Rubberized Surface For Comfortable Typing
- Supports Linux, Win 11, 10, 8, 7, Vista® & XP
- Rigid, Rugged & Washable. Easy To Clean Surface.
- Medical Certification
- Harsh Industrial, Military & Medical Applications.
- Complete White Colour
- Custom Layout Options & Keys Available
- USB Interface



## HU-KBD SPECIFICATIONS

| PANEL                     |                                                     |                                                                                |
|---------------------------|-----------------------------------------------------|--------------------------------------------------------------------------------|
| Features                  | Keyboard                                            | Industrial Wireless Silicone Compact Keyboard With Function Keys.              |
|                           | Touchpad                                            | High dpi Work Area 55 x 40 mm Scroll Function - Left & Right Mouse Buttons.    |
|                           | Protection                                          | Ip67 Waterproof & Dustproof Protection Level.                                  |
|                           | Washable                                            | Easy To Clean Surface - Any Hospital-Approved Disinfectants.                   |
|                           | Application                                         | Harsh Industrial, Military, Medical Environments Where Harmful Bacteria Exists |
| Material:                 | Silicone                                            | Industrial, Military-Grade & Medical Grade Coating.                            |
|                           | Colour                                              | White                                                                          |
| Electrical Data           | Interface                                           | 2 x AAA Alkaline Batteries, Working: 15mA, Stand-By: 5mA, Sleep: 13uA          |
|                           | Working Distance                                    | -20° C to + 80° C, -30° C to + 80° C                                           |
|                           | Powered By, Power Consumption                       | Optional Antibacterial Component In Silicone Available On Request              |
| Mechanical                | Total Keys                                          | 104                                                                            |
|                           | Function Keys                                       | 12                                                                             |
|                           | Key Lifespan                                        | > 10 Million Actuations                                                        |
|                           | Language & Layout                                   | Standard US English (All Windows Compatible Languages Available On Request)    |
|                           | Other                                               | Custom Keys & Layouts Available                                                |
| Dimension                 | Product, Packaging                                  | 418(L) x 133(W) x 16(H) mm, 440(L) x 175(W) x 22(H) mm                         |
| Weight                    | Net, Gross                                          | 830g, 980g                                                                     |
| Operation Systems Support | Linux, Win 11, 10, 8, 7, Vista® & XP, Mac OS9 / OSX |                                                                                |