


OUR TRUSTED LEGACY IS YOUR SUCCESS



PRODUCT CATEGORY:
BEDHEAD SERVICE SYSTEMS

PRODUCT NAME:
HU20-01

PRODUCT DESCRIPTION:
VERTICAL BEDHEAD SERVICE SYSTEM



Area of Use: Normal Care, Neonatal, Maternity, Medical Consultation Rooms, Recovery Wards and Trauma Units

Intended Use: The bedhead service unit is fixed next to the patient bed or trolley position, and is intended to supply one or more of the following services:

- electrical power up to 240 volts
- medical gas supply - oxygen, nitrous oxide, carbon dioxide, vacuum
- low pressure medical air for patient ventilation
- high pressure medical air for driving other medical devices
- lighting for illumination purposes
- communication services

for diagnosis and communication purposes in medical areas and special areas. The equipment is intended to be used by qualified nursing staff, Medical Practitioners and patients (e.g. nurse call, light).

Beds: Multiple bed application.

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU20-01 SPECIFICATIONS

SPECIFICATIONS:	
Construction:	Aluminium extrusion.
Surface Finish:	Natural Anodised 25 μ , powder coating in various colours or vinyl wrapping for centre fascia panel.
Installation Position:	Vertical Wall Mounted next to the patient's bed or sharing between beds.
Service Entry:	Via the ceiling void.
High Voltage:	Various switch sockets can be provided to meet the specific requirements.
Low Voltage:	Installed to suit specific requirements, e.g.: data, telephone, network plugs, earth leakage.
Medical Gas:	Various types can be installed or provision made for own installation.
Illumination:	1425 lumen LED Reading Light (15W, 4000K), 950 lumen LED Nightlight (10W, 4000K)
Customised Option:	Unit length and service requirements are manufactured to meet the specific needs.
Product Weight:	15 Kg/m
Payload of Equipment Pole:	20 kg per side (These payloads pertain to brick/concrete walls. Drywall payloads to be supplied on demand.)
Suitable Accessories:	Medical Rails, Infusion Bag Holders, Examination Light, Monitor Support Systems, Equipment Trays and Equipment Clamps.
Certification:	ISO13485



TECHNICAL DRAWINGS

