

ADVANCING HEALTHCARE TECHNOLOGY



PRODUCT CATEGORY:

MEDICAL ACCESSORIES - ATTACHMENTS

PRODUCT NAME:

HU-PET

PRODUCT DESCRIPTION:

PLASTIC EQUIPMENT TRAY



Area of Use: All Areas

Intended Use: The HU-PET range of plastic equipment trays were specifically made to provide more space to store medical consumables, instruments or other medical equipment in. The unit comes in various mounting options to attach to our Hutz Medical range of Pendant Service Systems and Bedhead Service Systems.

ADVANCING HEALTHCARE TECHNOLOGY



HU-PET SPECIFICATIONS

SPECIFICATIONS:	
Description:	The HU-PET is a versatile storage container with a solid powder coated frame. The container can be lifted out of the frame for easy cleaning.
Construction:	Aluminium
Surface Finish:	Power coated (RAL 9002) as a standard.
Payload:	10 Kg



HU-PET CONFIGURATIONS:

PRODUCT CODE:	PRODUCT IMAGE:	TECHNICAL DRAWING	INSTALLATION POSITION	DIM X (MM)
HU-PET-P-01			Pole Mounted	321 mm x 161 mm x 148 mm
COMPATABILITY: All Hutz Medical Pendant and Bedhead Service System.				
HU-PET-W-01			Wall Mounted	321 mm x 161 mm x 148 mm
HU-PET-R-01			Rail Mounted	321 mm x 161 mm x 148 mm
COMPATABILITY: All Hutz Medical Pendant and Bedhead Service Systems / Gabler Rails.				