

ADVANCING HEALTHCARE TECHNOLOGY



PRODUCT CATEGORY:

MEDICAL ACCESSORIES - ATTACHMENTS

PRODUCT NAME:

HU-WEB

PRODUCT DESCRIPTION

WIRE EQUIPMENT BASKET



Area of Use: All Areas

Intended Use: Wire Equipment baskets were specifically made to provide more space to store medical consumables, instruments or other medical equipment in.

ADVANCING HEALTHCARE TECHNOLOGY



HU-WEB SPECIFICATIONS

SPECIFICATIONS:	
Construction:	Stainless Steel
Surface Finish:	Electro-polished
Installation Position:	Wall / Pole / Rail Mounted
Grid Size:	Ø3 mm @ 20 mm x 30 mm pitch / Ø3 mm @ 25 mm x 55 mm pitch



HU-WEB CONFIGURATIONS:

PRODUCT CODE:	PRODUCT IMAGE:	TECHNICAL DRAWING	CONFIGURATION	CONSTRUCTION	DIM x (mm)
HU-WEB-01			Wall / Pole / Rail Mounted With or without handle	Stainless Steel	Various