

ADVANCING HEALTHCARE TECHNOLOGY



PRODUCT CATEGORY:

PENDANT SERVICE SYSTEMS

PRODUCT NAME:

HU8-FSP-01

PRODUCT DESCRIPTION:

FLOOR STANDING PEDESTAL



Area of Use: Casualty, Trauma and Emergency Units.

Intended Use: The pendant service system is floor mounted and is intended to supply one or more of the following medical services via a fixed, service pillar:

- electrical power up to 240 volts
- medical gases such as oxygen, nitrous oxide, carbon dioxide
- low pressure medical air for patient ventilation
- high pressure medical air for driving other medical devices
- vacuum and / or anaesthetic gas scavenging services

The unit can be fitted with medical rails for monitors, life support systems, medical attachments etc.

The equipment is intended to be used by surgeons, anaesthetists, qualified Trauma, Theatre and Intensive / High Care nursing staff.



HU8-FSP-01 SPECIFICATIONS

SPECIFICATIONS:	
Construction:	The compact service column is manufactured from aluminium extrusion.
Surface Finish:	As a standard, fascia plates are offered in hard wearing textured powder coated RAL9002 finish or vinyl wrapping.
Installation Position:	Vertically mounted, floor fixed.
Service Entry:	Via a core drilled hole or chased in under the floor channels with services supplied from the plant room (Services enter through the base of the unit).
High Voltage:	Various switch sockets can be provided to meet the specific requirements.
Low Voltage:	Installed to suit specific requirements, e.g.: data, telephone, network plugs, earth leakage.
Medical Gas:	Various types can be installed or provision made for own installation.
Illumination:	Not applicable to this unit.
Customised Option:	Unit length and service requirement are manufactured to meet specific needs.
Rotation:	320° (Table Top Only)
Product Payload:	20 Kg
Product Weight:	45 Kg
Suitable Accessories:	Medical Rail, Drip-hangers, Plastic Equipment Tray, Wire Basket, Examination Light, Monitor Support System and Equipment Clamps.



TECHNICAL DRAWINGS

