

ADVANCING HEALTHCARE TECHNOLOGY



PRODUCT CATEGORY:
MEDICAL ACCESSORIES

PRODUCT NAME:
HU-IBH-M

PRODUCT DESCRIPTION:
MOBILE INFUSION BAG HOLDERS



Area of Use: All Areas

Intended Use: The HU-IBH Infusion Bag Holders are used to hold bags of fluid including medicine that are being administered to the patient at a steady flow. This is done by hanging a bag of fluid on the hook near the top of the pole. The unit is intended to be used by Medical Practitioners and Nursing Staff.



HU-IBH-M SPECIFICATIONS

SPECIFICATIONS:	
Description:	The HU-IBH Infusion Bag Holders are designed to hold bags of fluid that are being administered to the patient at a steady flow. They are made from stainless steel/aluminium with stainless steel hooks in various configurations.
Installation position:	Mobile




HU-IBH-M CONFIGURATIONS:

PRODUCT CODE:	PRODUCT IMAGE:	CONFIGURATION	CONSTRUCTION	HOOK ASSEMBLY	BASE	CASTORS
HU-IBH-M-01		4 Hooks With stainless steel screw in bottle holder. One hand height adjustment.	Stainless Steel	Stainless Steel	5 Star	5 x 50 mm twin castors, two with locks.
HU-IBH-M-02		4 Hooks Height adjustment via screw mechanism.	Stainless Steel	Stainless Steel	5 Star	5 x 60mm twin castors, two with locks.
HU-IBH-M-06		4 Hooks One hand height adjustment.	Robust Chrome	ABS Plastic	5 Star	5 x 50mm twin castors, two with locks.



HU-IBH-M CONFIGURATIONS:

PRODUCT CODE:	PRODUCT IMAGE:	CONFIGURATION	CONSTRUCTION	HOOK ASSEMBLY	BASE	CASTORS
HU-IBH-M-07		4 Hooks One hand height adjustment.	Stainless Steel	Stainless Steel	5 Star	5 x 80mm castors, two with locks.
HU-IBH-M-08		4 Hooks One hand height adjustment.	Robust Chrome	ABS Plastic	5 Star	5 x character based castors, two with locks.
HU-IBH-M-15		2 Hooks Height adjustment via screw mechanism (max height 1320mm / min height 830mm).	Black Coated Steel	Stainless Steel	3 Leg Stand Ø660 mm	Foldable legs on stand base.