

ADVANCING HEALTHCARE TECHNOLOGY



PRODUCT CATEGORY:

MEDICAL ACCESSORIES - ATTACHMENTS

PRODUCT NAME:

HU-IBH-P

PRODUCT DESCRIPTION:

POLE MOUNTED INFUSION BAG HOLDERS



Area of Use: All Areas

Intended Use: The HU-IBH Infusion Bag Holders are used to hold bags of fluid including medicine that are being administered to the patient at a steady flow. This is done by hanging a bag of fluid on the hook near the top of the pole. The unit is intended to be used by Medical Practitioners and Nursing staff.



HU-IBH-P SPECIFICATIONS

SPECIFICATIONS:	
Description:	The HU-IBH Infusion Bag Holders are designed to hold bags of fluid that are being administrated to the patient at a steady flow. They are made from stainless steel with stainless steel hooks in various configurations.
Installation position:	Pole Mounted



HU-IBH-P CONFIGURATIONS:

PRODUCT CODE:	PRODUCT IMAGE:	TECHNICAL DRAWING	CON-FIGURATION	CON-STRUCTION	HOOK ASSEMBLY	POLE CLAMP	POLE LENGTH
HU-IBH-P-01			2 Hooks Vertical T-Shape	Stainless Steel Ø16 mm pole	Stainless Steel 304	Ø38 mm clamp Aluminium	500 mm
HU-IBH-P-02			4 Hooks Vertical T-Shape	Stainless Steel Ø16 mm pole	Stainless Steel 304	Ø38 mm clamp Aluminium	500 mm
HU-IBH-P-03			6 Hooks Vertical T-Shape	Stainless Steel Ø16 mm pole	Stainless Steel 304	Ø38 mm clamp Aluminium	500 mm
HU-IBH-P-04			4 Hooks Horizontal	Stainless Steel Ø16 mm pole	Stainless Steel 304	Ø38 mm clamp Aluminium	306.5 mm
HU-IBH-P-05			8 Hooks Horizontal	Stainless Steel Ø16 mm pole	Stainless Steel 304	Ø38 mm clamp Aluminium	500 mm