



PRODUCT CATEGORY:

MEDICAL ACCESSORIES - ATTACHMENTS

PRODUCT NAME:

HU-COS

PRODUCT DESCRIPTION:

POLE MOUNTED CABLE ORGANISER



Area of Use: All Areas

Intended Use: The HU-COS is a cable organiser for the neat arrangement of cables and IV lines on the Hutz Medical range of Pendant and Bedhead Service Systems.

ADVANCING HEALTHCARE TECHNOLOGY




HU-COS SPECIFICATIONS

SPECIFICATIONS:	
Description:	The HU-COS is a Pole mounted injection moulded cable clip designed to fit the standard Ø38mm equipment pole.
Installation position:	Pole Mounted.
Configuration:	HU-COS-P Cable organisers for Ø38mm equipment poles.



HU-COS CONFIGURATIONS:

PRODUCT CODE:	PRODUCT IMAGE:	TECHNICAL DRAWING	CONFIGURATION	DIM X (MM)
HU-COS-P			Cable organiser for Ø 38 mm equipment poles.	450 mm x 448 mm
COMPATABILITY: All Hutz Medical Pendant and Bedhead Service Systems.				