

ADVANCING HEALTHCARE TECHNOLOGY



PRODUCT CATEGORY:  
**MEDICAL ACCESSORIES**

PRODUCT NAME:  
**HU-MER**

PRODUCT DESCRIPTION:  
**MEDICAL EQUIPMENT RAILS**



**Area of Use:** All Areas

**Intended Use:** The HU-MER is a medical equipment rail that provides a universal mounting system for the temporary mounting and support of any patient care equipment required within a healthcare environment.

ADVANCING HEALTHCARE TECHNOLOGY



## HU-MER SPECIFICATIONS

| SPECIFICATIONS:        |                                                                                                          |
|------------------------|----------------------------------------------------------------------------------------------------------|
| Description:           | The HU-MER is a medical equipment rail that is either wall mounted or to various Hutz Medical Equipment. |
| Construction:          | Anodised Aluminium.                                                                                      |
| Installation position: | Wall / Pole / Unit Mounted.                                                                              |
| Mounting Accessories:  | Supplied with pole mounted clamp.                                                                        |
| Payload:               | > 20 Kg                                                                                                  |
| Shape:                 | Gabler shape / Euro shape.                                                                               |



## HU-MER CONFIGURATIONS:

| PRODUCT CODE:       | PRODUCT IMAGE: | TECHNICAL DRAWING | CON-FIGURATION                                      | DIM X (MM)      | PRODUCT COMPATABILITY                                  |
|---------------------|----------------|-------------------|-----------------------------------------------------|-----------------|--------------------------------------------------------|
| POLE MOUNTED        |                |                   |                                                     |                 |                                                        |
| HU-MER-P-01/02      |                |                   | Gabler Shape or Euro Shape<br><br>(Ø25 mm / Ø38 mm) | 150 mm          | All Hutz Medical Pendants and Bedhead Service Systems. |
| HU-MER-P-03/04      |                |                   |                                                     | 300 mm          |                                                        |
| HU-MER-P-05/06      |                |                   |                                                     | 400 mm          |                                                        |
| HU-MER-P-07/08      |                |                   |                                                     | 600 mm          |                                                        |
| DOUBLE POLE MOUNTED |                |                   |                                                     |                 |                                                        |
| HU-MER-DP-01        |                |                   | Gabler Shape or Euro Shape<br><br>(Ø25 mm / Ø38 mm) | Various Lengths | All Hutz Medical Pendants and Systems.                 |
| HU-MER-DP-02        |                |                   |                                                     |                 |                                                        |
| HU-MER-DP-475-01    |                |                   |                                                     |                 |                                                        |
| HU-MER-DP-475-02    |                |                   |                                                     |                 |                                                        |
| WALL MOUNTED        |                |                   |                                                     |                 |                                                        |
| HU-MER-W-01         |                |                   | Gabler Shape or Euro Shape                          | Various Lengths | Wall                                                   |
| DRY WALL MOUNTED    |                |                   |                                                     |                 |                                                        |
| HU-MER-DW-01        |                |                   | Gabler Shape or Euro Shape                          | Various Lengths | Dry Wall                                               |
| HU-MER-DW-02        |                |                   |                                                     |                 |                                                        |
| UNIT MOUNTED        |                |                   |                                                     |                 |                                                        |
| HU-MER-U-01         |                |                   | Gabler Shape or Euro Shape                          | Various Lengths | All Hutz Medical Bedhead Service Systems.              |
| HU-MER-U-02         |                |                   |                                                     |                 |                                                        |