

ADVANCING HEALTHCARE TECHNOLOGY



PRODUCT CATEGORY:
MEDICAL ACCESSORIES

PRODUCT NAME:
HU-IBH-R

PRODUCT DESCRIPTION:
RAIL MOUNTED INFUSION BAG HOLDERS



Area of Use: All Areas

Intended Use: The HU-IBH Infusion Bag Holders are used to hold bags of fluid including medicine that are being administered to the patient at a steady flow. This is done by hanging a bag of fluid on the hook near the top of the pole. The unit is intended to be used by Medical Practitioners and Nursing staff.



HU-IBH-R SPECIFICATIONS

SPECIFICATIONS:	
Description:	The HU-IBH Infusion Bag Holders are designed to hold bags of fluid that are being administrated to the patient at a steady flow. They are made from stainless steel/aluminium with stainless steel hooks in various configurations.
Installation position:	Rail Mounted



HU-IBH-R CONFIGURATIONS:

PRODUCT CODE:	PRODUCT IMAGE:	TECHNICAL DRAWING	CON-FIGURATION	CON-STRUCTION	HOOK ASSEMBLY	POLE CLAMP	POLE LENGTH
HU-IBH-R-09			2 Hooks Fixed Unit Vertical T-Shape	Aluminium Ø25 mm pole	Stainless Steel 304	Ø25 mm clamp Aluminium	295 mm
HU-IBH-R-10			4 Hooks Fixed Unit Vertical T-Shape	Aluminium Ø25 mm pole	Stainless Steel 304	Ø25 mm clamp Aluminium	375 mm
HU-IBH-R-11			3 Hooks Adjustable Unit Vertical T-Shape	Stainless Steel Ø25 mm pole	Stainless Steel 304	Ø40 mm clamp Aluminium	1150 mm - 1700 mm
HU-IBH-R-12			4 Hooks Adjustable Unit Vertical T-Shape	Stainless Steel Ø25 mm pole	Stainless Steel 304	Ø40 mm clamp Aluminium	1150 mm - 1700 mm

ADVANCING HEALTHCARE TECHNOLOGY



HU-IBH-R CONFIGURATIONS:

HU-IBH-R-13			3 Hooks Adjustable Unit Angled with Wall Stay	Aluminium extruded Ø16mm pole	Stainless Steel 304	Ø40 mm clamp Aluminium	1180 mm - 1501 mm
HU-IBH-R-16			3 Hooks Fixed Unit Vertical T-Shape	Aluminium extruded Ø16mm pole	Stainless Steel 304	Ø40 mm clamp Aluminium	625 mm