

ADVANCING HEALTHCARE TECHNOLOGY



PRODUCT CATEGORY:
MEDICAL ACCESSORIES

PRODUCT NAME:
HU-MEP

PRODUCT DESCRIPTION:
MEDICAL EQUIPMENT POLES



Area of Use: All Areas

Intended Use: The HU-MEP is a medical equipment pole that provides a universal mounting system for the temporary mounting and support of any patient care equipment required within a healthcare environment.

ADVANCING HEALTHCARE TECHNOLOGY



HU-MEP SPECIFICATIONS

SPECIFICATIONS:	
Description:	The HU-MEP is a medical equipment pole that provides a universal mounting system for the temporary mounting and support of any patient care equipment.
Construction:	304 Stainless Steel Pole.




HU-MEP CONFIGURATIONS

PRODUCT CODE:	PRODUCT IMAGE:	CON-FIGURATION	PAYLOAD:	DIAMETER:	DIM X (MM)	MOUNTING ACCESSORIES:	PRODUCT COMPATABILITY
HU-MEP-W		Wall Mounted	> 20 Kg	Ø 25 mm or Ø 38 mm	Various lengths	Supplied with spacers & anchors.	N/A
HU-MEP-DW		Dry Wall Mounted	> 20 Kg	Ø 38 mm	1000 mm	Supplied with spacers & anchors.	N/A
HU-MEP-R-01		Rail Mounted	> 20 Kg	Ø25mm or Ø 38mm	1000 mm	Supplied with single clamp and stay.	All Hutz Medical Bedhead Service Systems / Medical Rails.
HU-MEP-R-03					800 mm		
HU-MEP-R-02		Rail Mounted	> 20 Kg	Ø25mm or Ø 38mm	1000 mm	Supplied with dual clamps. Supplied with 2 x 75mm clamps.	All Hutz Medical Bedhead Service Systems / Medical Rails.
HU-MEP-R-05				Ø 38mm			



HU-MEP CONFIGURATIONS

PRODUCT CODE:	PRODUCT IMAGE:	CON-FIGURATION	PAYLOAD:	DIAMETER:	DIM X (MM)	MOUNTING ACCESSORIES:	PRODUCT COMPATABILITY
HU-MEP-E-01		Pole Mounted	> 20 Kg	Ø 25 mm clamp to Ø 25mm pole.	1000 mm	Supplied with 2 x pole holder extender clamps made from anodised aluminium.	All Hutz Medical Pendant Service Systems and HU20 Bedhead Service Systems.
HU-MEP-E-02				Ø 25 mm clamp to Ø 38mm pole.			
HU-MEP-E-04		Pole Mounted	> 20 Kg	Ø 25 mm clamp to Ø 25mm pole.	1000 mm	Supplied with 2 x pole holder extender clamps made from anodised aluminium.	All Hutz Medical Pendant Service Systems and HU20 Bedhead Service Systems.