

ADVANCING HEALTHCARE TECHNOLOGY



PRODUCT CATEGORY:

MEDICAL ACCESSORIES - ATTACHMENTS

PRODUCT NAME:

HU-EH

PRODUCT DESCRIPTION:

EQUIPMENT HOOK



Area of Use: All Areas

Intended Use: The HU-EH is a small hook made from stainless steel, which serves to grab, latch or in any way attach itself onto another object, used to store/hang equipment like scopes etc. The unit comes in various mounting options to attach to our Hutz Range of Pendant Service Systems and Bedhead Service Systems.

ADVANCING HEALTHCARE TECHNOLOGY



HU-EH SPECIFICATIONS

SPECIFICATIONS:	
Description:	The HU-EH is a small hook made from stainless steel, which serves to grab, latch or in any way attach itself onto another object, used to store/hang equipment like scopes etc. The unit comes in various mounting options to attach to our Hutz Range of Pendant Service Systems and Bedhead Service Systems.
Installation position:	Pole and Rail Mounted



HU-ABH CONFIGURATIONS:

PRODUCT CODE:	PRODUCT IMAGE:	TECHNICAL DRAWING	CON-FIGURATION	CON-STRUCTION	DIM x (mm)	Clamp
HU-EH-P			Pole Mounted	Stainless Steel	181.9mm x 100mm x 72.5mm	Supplied with an aluminium clamp to suit a Ø25 mm and Ø38 mm equipment pole.

COMPATABILITY: All Hutz Medical Pendants and Bedhead Service Systems.

HU-EH-R			Rail Mounted	Stainless Steel	171.8mm x 85.0mm x 65.5mm	Supplied with an aluminium clamp to suit a Gabler medical rail.
---------	--	--	--------------	-----------------	---------------------------	---

COMPATABILITY: All Hutz Medical Pendants and Bedhead Service Systems.