

OUR TRUSTED LEGACY IS YOUR SUCCESS



PRODUCT CATEGORY:
MEDICAL ACCESSORIES

PRODUCT NAME:
HU-T-MA-M-03

PRODUCT DESCRIPTION:
**DUAL CEILING MOUNTED MONITOR ARM FOR 4 MONITORS
WITH BEZEL PROTECTION**



Area of Use: All Areas

Intended Use: The Arms for Monitors are ceiling mounted systems specially designed for operating theatre. The monitors can be well fixed on the extension flexible spring arm which can be easily moved, rotated and height adjusted to the optimum position and ergonomically facilitate working procedures for operation.

OUR TRUSTED LEGACY IS YOUR SUCCESS

HU-T-MA-M-03 SPECIFICATIONS

SPECIFICATIONS:	
Construction:	The arm and pillar are manufactured from aluminium extrusion. Main body is made by high strength extrusion aluminum alloy and steel for lightweight and robust.
Surface Finish:	Surface with mold plastic cover, anti-corrosion, easy to clean and disinfect.
Installation Position:	Ceiling suspended.
Service Entry:	Via a ceiling void or a core hole with services supplied from the plant room.
Vertical Lift:	≤600 mm
Motor Power:	≤1KW
Brake System:	Dual brake system of mechanical friction brake and pneumatic brake with control can keep the pendant from drifting and prevent unintentional movement which may cause disconnection of life supporting systems.
Rotation:	Swiveling arm with 340 ° degree rotation.
Product Payload:	Up to 100 Kg
Load Rating:	400 Kg
Rated Voltage:	220V; Rated frequency: 50HZ; Input power: 4KVA;
Console:	Console with monitor fixing bracket with bezel protection(Monitor is not included) x 1 set
High Voltage:	3 pcs x power sockets
	2 pcs x BNC socket
Frame:	Bezel frame with control x 1 pc
Data:	Fiber optic cable, RJ45 internet cable, VGA cable should be provided by monitor supplier which should thread to pendant console and connect with monitor directly.
Certification:	CE

ADDITIONAL INFORMATION:

