



PRODUCT CATEGORY:
HOSPITAL FURNITURE

PRODUCT NAME:
HU-U0035.YL3

PRODUCT DESCRIPTION:
HOSPITAL OVERBED TABLE



Product Codes: HU-U0035.YL3

Area of Use: Normal Care to ICU Wards

Intended Use: An overbed table is utilised by both patients and healthcare workers alike. They are a surface on which to place food trays and can hold the patient's personal items such as books or phones; while healthcare workers utilise them to hold medical devices or supplies.




HU-U0035.YL3 SPECIFICATIONS

SPECIFICATIONS:	
Table Top:	Compact Laminate top that is durable against water and heat.
Height Adjustment:	Height of overbed table is adjusted by gas spring according to height of bed.
Material:	Epoxy painted steel construction, resistant to corrosion.
Material - Leg:	4 x castors, two with brakes.
Width:	800 mm
Depth:	400 mm
Height: (Adjustable)	790 mm - 1030 mm
Weight:	+/-10,5 kg



LAMINATE COLOUR OPTIONS & CODES

 RL01 : GENTAS 3156	 RL02 : GENTAS 3132	 RL04 : GENTAS 3104	 GENTAS 3153	 RL06 : GENTAS 3155
 RL07 : GENTAS 3120	 RL08 : GENTAS 3134	 RL10 : GENTAS 3146	 RL11 : GENTAS 3147	 RL12 : GENTAS 3151
 RL14 : GENTAS 4206	 RL15 : GENTAS 4251	 RL16 : GENTAS 4212	 RL17 : GENTAS 4234	