



PRODUCT CATEGORY:
HOSPITAL FURNITURE

PRODUCT NAME:
HU-I-U0041

PRODUCT DESCRIPTION:
EMERGENCY STRETCHER



Area of Use: Resus Area / Trauma and Emergency

Intended Use: A stretcher is an apparatus used for moving patients who require medical care. A wheeled stretcher is often equipped with variable height frames, wheels, tracks, or skids. Stretchers are primarily used in acute out-of-hospital care situations by Emergency Medical Services (EMS), military, and Search and Rescue personnel.





HU-I-U0041 SPECIFICATIONS

SPECIFICATIONS:	
Area:	Resus Area / Trauma and Emergency
Base Material:	Easy cleanable lower frame covered by ABS Plastic.
Metal Components:	All metal components are powder coated.
Mattress Platform:	Easy to clean, x-ray translucent compact laminate (HPL) mattress platform.
Bumpers:	There are protective bumpers on all four corners of the bed.
Trendelenburg / Reverse Trendelenburg:	0 - 13° / 0 - 13° (Hydraulic pedal system can be controlled from both sides).
Adjustable Backrest:	0 - 62°
Height Adjustment:	690 mm - 1010 mm (Dual sided foot controlled hydraulic system for raising and lowering the bed).
Castors:	4 x 150 mm diameter castors / Central locking brake system that can be (de)activated from either end of the bed.
Rails:	Adjustable side rails that release at the base of the bed for maximum safety.
Dimensions:	2230 mm (Length) / 780 mm (Width)
Safe Working Load:	310 Kg
Cleaning:	Power washable bed and mattress.
Mattress Thickness:	Two piece, 80 mm thick, waterproof mattress.
Mattress Density:	28 Kg/m ³
Mattress Cover:	Water resistant, antibacterial, breathable imitation leather cover.
IV Pole Provision:	The stretcher makes allowance for the use of IV poles on all four corners.
Volume:	1,3 m ³
Accessories supplied with:	Integrated Oxygen Cylinder Holder.
	Storage Tray on the base.
Options:	Vinlex colour choice.



OPTION CODES

GROUP 1 / CASTOR OPTIONS



GROUP 2 / IV STAND OPTIONS



DESCRIPTION OF CODES:

- T11** CASTOR TWO WITH BRAKE (150 MM)
- T12** CASTOR CENTRAL & DIRECTIONAL LOCK (150 MM)
- T17** CASTOR CENTRAL & DIRECTIONAL LOCK (150 MM)
- W0** WITHOUT OTHER ACCESSORIES
- W8** IV STAND : ONE HEIGHT ADJUSTMENT SYSTEM

Must select only one of them.



COLOUR OPTION CODES

