



PRODUCT CATEGORY:
HOSPITAL FURNITURE

PRODUCT NAME:
HU-I-U0138

PRODUCT DESCRIPTION:
ICU CHART TROLLEY



Product Codes: HU-I-U0138

Area of Use: ICU and CCU

Intended Use: A Chart Trolley is utilised for the storage of medication, supplies and to write patients' charts at the end of the bed in an ICU or CCU. The Chart Trolley also includes a patient storage cabinet purposely built to stand at the end of every bed.



HU-I-U0138 SPECIFICATIONS

SPECIFICATIONS:	
Description:	Chart Trolley for storage of medication, supplies and to write patients' charts at the end of the bed in an ICU or CCU. It also includes a patient storage cabinet purposely built to stand at the end of every bed. The unit features a smooth, refined and hardwearing design.
Trolley Material:	Manufactured from steel with a white epoxy coated finish.
Trolley Frame Material:	The frame is manufactured from 30 x 30 mm square tubing.
Trolley Top Material:	The sloping top is manufactured from covered laminated chipboard allowing a comfortable writing surface for a large ICU chart, while still ensuring visibility of the patient.
Trolley Top Lip:	32 mm Aluminum lip at the bottom of writing surface to prevent the chart from slipping.
Dimensions:	
Length - Top:	1150 mm
Width - Top:	600 mm
Height - Floor to Top:	1000 mm
Weight:	40 kg
Shelf:	1 x Large shelf underneath the writing surface (top), closed on three sides.
Drawers:	4 x Pull out drawers on the left-hand side of the unit.
Patient Storage:	1 x Pull out drawer as well as 1 x locker with self inside and door catch on the right-hand side.
Castors:	4 x Grey rubber 75 mm castors (2 x lockable and 2 x swivel).
File/Chart Holder:	The unit is supplied with a file and X-ray chart board holder on left-hand side.