



PRODUCT CATEGORY:
HOSPITAL FURNITURE

PRODUCT NAME:
HU-R-01

PRODUCT DESCRIPTION:
RECLINER CHAIR



Area of Use: Maternity Ward / Paediatrics / Oncology / General Wards / Medical Wards /Surgical Wards

Intended Use: Reclining chairs are used in medical facilities to keep patients comfortable while recovering from illness or surgery or undergoing procedures. They are built to be durable, functional and low maintenance. The dedicated chairs are specially designed to support and add comfort , with lumbar support, armrests and reclining features.



HU-R-01 SPECIFICATION

SPECIFICATIONS:	
Description:	Comfortable, durable recliner for nursing infants. With enough support to allow the mother to stand up easily after nursing without disturbing the infant.
Upholstery:	Synthetic leather.
Positions:	Two positions: Seated and Reclined.
Additional Support:	Adjustable Headrest, armrests and footrest.
Base Material:	Epoxy steel.
Feet:	Black Plastic Stumps with protective ends.
Width (External):	700 mm
Width (Sitting Area):	535 mm
Depth:	950 mm – 1800 mm (upright to reclined)
Height (Sitting):	460 mm
Height (Backrest):	1200 mm
Weight:	120 Kg
Optional:	Castors with Brakes.
	Side Pockets