

ADVANCING HEALTHCARE TECHNOLOGY



PRODUCT CATEGORY:

MEDICAL AND SURGICAL LIGHTING

PRODUCT NAME:

LUVIS E100

PRODUCT DESCRIPTION:

MEDICAL EXAMINATION/PROCEDURE LIGHT



Product Codes: LUVIS E100

Area of Use: Normal Care to High Care Wards, Medical Consultation Rooms, Intensive Care Units, Trauma Units, Nursing Homes and Clinics.

Intended Use: The Luvivis E100 Medical Examination/Procedure Light is designed to provide illumination of the body of a patient for localized diagnosis and treatment. The light is not designed for use during surgery. It is designed for continuous operation and it must be possible to stop the treatment or diagnosis without danger to the patient if the light fails during use. The device is intended to be used by doctors / medical practitioners and nursing staff.

ADVANCING HEALTHCARE TECHNOLOGY



LUVIS E100 SPECIFICATIONS

SPECIFICATIONS:	
Multicoloured Light:	Yes
Central illuminance at 500 mm:	60 000 Lux @ 500 mm
Colour Temperature:	3500 K / 4 000 K / 4500 K - 3 steps
Colour rendering index (Ra):	95
Number of LED's:	8
LED Lifespan:	>50 000 Hrs
Head size:	250 mm x 246 mm
Mounting Options:	Wall /Mobile /Ceiling Mounted
Power Consumption:	17W
Intensity Control:	3 Steps
Certification:	CE
Additional Options:	
Wall Control Option:	N/A
Battery Backup Option:	Yes (Mobile)



OTHER PICTURES:



Images are used for illustration purposes only.