

OUR TRUSTED LEGACY IS YOUR SUCCESS



PRODUCT CATEGORY:
VETERINARY LIGHTING

PRODUCT NAME:
MASTIFF

PRODUCT DESCRIPTION:
VETERINARY PROCEDURE LIGHT



Area of Use: Veterinary Rooms / Veterinary Theatres / Veterinary Clinics

Intended Use: The Mastiff Veterinary Procedure Light is designed to provide illumination for localised diagnosis and treatment. The light is versatile featuring a three-lighting-mode function. The device is intended to be used by Veterinarians / Veterinarian staff.

OUR TRUSTED LEGACY IS YOUR SUCCESS



MASTIFF SPECIFICATIONS

SPECIFICATIONS:	
Two-Lighting-Mode Function:	Surgical (Default) / Examination-Procedure (Activated with button)
Central illuminance:	130 000 lux (Examination-Procedure: 60 000 lux)
Intensity Control from/to:	40 000 - 130 000 lux - 3 steps
Adjustable Focus:	150 - 250 mm (Surgical Mode) / Examination-Procedure Mode: 150 mm
Colour Temperature:	3800 K / 4 300 K / 4800 K - 3 steps (Dental Light Mode: 4800 K)
Colour rendering index (Ra/R9):	96 / 95
Number of LED's:	42
LED Lifespan:	>50 000 Hrs
Head size:	435 mm
Mounting Options:	Ceiling Mounted (Mastiff-C) / Wall Mounted (Mastiff-W)
Detachable Centre Grip:	Yes (x2)
Certification:	CE
Additional Options:	
Wall Control Option:	Yes (Tablet PC)
Battery Backup Option:	Yes
Weight:	Single: 28 Kg



OTHER CONFIGURATIONS:



Images are used for illustration purposes only.