

ADVANCING HEALTHCARE TECHNOLOGY



PRODUCT CATEGORY:
MEDICAL ACCESSORIES

PRODUCT NAME:
HU-MSS

PRODUCT DESCRIPTION:
MONITOR SUPPORT SYSTEMS



Area of Use: All Areas

Intended Use: The HU-MSS range of monitor brackets allows the nursing staff to mount a computer screen or patient monitor to larger medical equipment such as a pendant or bedhead service system.



HU-MSS SPECIFICATIONS

SPECIFICATIONS:	
Description:	The HU-MSS is a universal monitor bracket with 110° adjustment to suit either base or rear mounting of the monitor. The unit also allows for tilt and swivel function to position the equipment.
Surface Finish:	Powder Coated
Payload:	14 -27 Kg (depending on front-end mounting configuration)
Installation Position:	Pole, Rail or Wall Mounted



HU-MSS CONFIGURATIONS:

PRODUCT CODE:	PRODUCT IMAGE:	CONFIGURATION	CONSTRUCTION
HU-MSS-W-01/02		Wall Mounted with or without basket	Casted Aluminium Arm
HU-MSS-P-01/02		Pole Mounted with or without basket	Casted Aluminium Arm
HU-MSS-P-03		Dual Pole Mounted with or without cable management with or without wire basket	Casted Aluminium Arm
HU-MSS-R-01/02		Rail Mounted with or without wire basket	Casted Aluminium Arm