

OUR TRUSTED LEGACY IS YOUR SUCCESS



PRODUCT CATEGORY:
MEDICAL ACCESSORIES

PRODUCT NAME:
HU-T-MA

PRODUCT DESCRIPTION:
CEILING/WALL MOUNTED MONITOR ARMS



Area of Use: All Areas

Intended Use: The Arms for Monitors are ceiling mounted systems specially designed for operating theatre. The monitors can be well fixed on the extension flexible spring arm which can be easily moved, rotated and height adjusted to the optimum position and ergonomically facilitate working procedures for operation.

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-T-MA SPECIFICATIONS

SPECIFICATIONS:	
Construction:	The arm and pillar are manufactured from aluminium extrusion. Main body is made by high strength extrusion aluminum alloy and steel for lightweight and robust.
Surface Finish:	Surface with mold plastic cover, anti-corrosion, easy to clean and disinfect.
Installation Position:	Ceiling suspended / Wall Mounted.
Service Entry:	Via a ceiling void or a core hole with services supplied from the plant room.
Rotation:	Swiveling arm with 340 ° degree rotation.
Rated Voltage:	220V; Rated frequency: 50HZ; Input power: 4KVA;
Control Handle:	1 pcs per arm.
Certification:	CE



HU-T-MA CONFIGURATIONS:

PRODUCT CODE:	PRODUCT IMAGE:	SCREEN CONFIGURATION	MOUNTING POSITION	LOAD CAPACITY	RADIUS / ADJUSTMENT	BRAKE	CONFIGURATION
HU-T-MA-S-01		Single Screen Monitor Arm	Ceiling Mounted	≤20 kg	Extension flexible spring arm with radius 850 mm + 910 mm Spring arm height adjustment range: 0-800mm	Mechanical friction brake	Monitor fixing bracket -1pc BNC point with 5m cable -1pc Power cable with 5m cable and plug -1pc
HU-T-MA-D-01		Double Screen Monitor Arm	Ceiling Mounted	≤20 kg	Extension flexible spring arm with radius 850 mm + 910 mm Spring arm height adjustment range: 0-800mm	Mechanical Friction brake	Monitor fixing bracket for double 27" monitor -1pc BNC point with 5 meter cable- 2pcs Power cable with 5 meter cable
HU-T-MA-D-02		Double Screen Double Arm	Ceiling Mounted	≤20 kg	2 x Extension flexible spring arm with radius 850 mm + 910 mm Spring arm height adjustment range: 0-800mm	Mechanical Friction brake	Each arm with a monitor fixing bracket for one monitor - 1pc. BNC point with 5m cable - 1pc for each arm Power cable with 5m cable and plug - 1pc for each arm

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-T-MA CONFIGURATIONS:

PRODUCT CODE:	PRODUCT IMAGE:	SCREEN CONFIGURATION	MOUNTING POSITION	LOAD CAPACITY	RADIUS / ADJUSTMENT	BRAKE	CONFIGURATION
HU-T-MA-D-03		Dual Screen Monitor Arm for 4 Monitors	Ceiling Mounted	≤20 kg	2 x extension flexible spring arm with radius 850mm + 910mm Spring arm height adjustment range: 0-800mm	Mechanical friction brake	Each arm with a monitor fixing bracket for one monitor - 1pc BNC point with 5m cable - 1pc for each arm Power cable with 5m cable - 1p for each arm
HU-T-MA-M-04		Single Screen Monitor Arm	Ceiling Mounted	≤20 kg	Double motorized arm with radius 750mm + 750mm Vertical adjustment/lift ≤430mm via ≤1KW motor power	Electro-pneumatic brake	Console with 55" large display monitor fixing bracket - 1pc Power socket - 3pcs (safely in concealed built in electrical box) Fibre optic cable, RJ45 internet cable, VGA cable should be provided by monitor supplier which should thread to pendant console and connect with monitor directly
HU-L-TA		Single Tablet Arm	Wall Mounted	≤4 kg	Full Extension 1137mm	-	Screen size 17"-32 Vesa compatible 75x75, 100x100