

OUR TRUSTED LEGACY IS YOUR SUCCESS



PRODUCT CATEGORY:

MEDICAL AND SURGICAL LIGHTING

PRODUCT NAME:

HUGO LED-D

PRODUCT DESCRIPTION:

MEDICAL EXAMINATION LIGHT



Area of Use:

Normal Care to High Care Wards, Medical Consultation Rooms, Intensive Care Units, Trauma Units, Nursing Homes and Clinics.

Intended Use:

The HUGO LED Medical Examination Light is designed to provide illumination of the body of a patient for localised diagnosis and treatment. The light is not designed for use during surgery. It is designed for continuous operation and it must be possible to stop the treatment or diagnosis without danger to the patient if the light fails during use. The device is intended to be used by doctors / medical practitioners and nursing staff.

OUR TRUSTED LEGACY IS YOUR SUCCESS




HUGO LED-D SPECIFICATIONS

SPECIFICATIONS:	HUGO LED-D (Desk Mounted)
Light Source:	5x LEDs
Luminance Output:	35 000 Lux @ 1000 mm / 50 000 Lux @ 500 mm
Supply Voltage:	220-240V/50-60Hz AC
Electrical Output:	12V DC/17.5 Watt
Electrical Safety Classification:	Class 2 Type B
Light Field:	150 mm @ 500mm / 180mm @ 1000mm
Colour Temperature:	4500 K
LED Lifespan:	>50 000 Hrs
Colour Rendering Index:	>85
Switching:	Illuminated switch on lamp head.
Lamp body/Housing: (All plastic Parts)	Plastic moulded.
Arm:	Aluminium: Natural anodised
Colour:	Lamp body - White, Arms - Natural anodised, cover caps offered in blue.
Mains Lead:	Moulded plug top, 1.5 m cable.
Product Weight (Net):	2.4 Kg
Dimensions (Packaging):	780 mm (H) x 200 mm (W) x 130 mm (D) *Folded
Operating Reach:	1200 mm
Mounting Options:	Desk / Table
Warranty:	2 Years
IP Rating:	54



TECHNICAL DRAWINGS



60000-70000
50000-60000
40000-50000
30000-40000
20000-30000
10000-20000
0-10000

