

OUR TRUSTED LEGACY IS YOUR SUCCESS



PRODUCT CATEGORY:
VETERINARY LIGHTING

PRODUCT NAME:
JACK RUSSELL

PRODUCT DESCRIPTION:
VETERINARY EXAMINATION LIGHT



Area of Use: Veterinary Rooms / Veterinary Theatres / Veterinary Clinics

Intended Use: The Jack Russell Veterinary Examination Light is designed to provide illumination of the body of a animal for localized diagnosis and treatment. The device is intended to be used by Veterinarians / Veterinarian staff.

OUR TRUSTED LEGACY IS YOUR SUCCESS



JACK RUSSELL SPECIFICATIONS

SPECIFICATIONS:	
Light Source:	5x LEDs
Luminance Output:	35 000 Lux @ 1000 mm / 50 000 Lux @ 500 mm
Supply Voltage:	220-240V/50-60Hz AC
Electrical Output:	12V DC/17.5 Watt
Electrical Safety Classification:	Class 2 Type B
Light Field @ 1000 mm:	Ø 180 mm
Colour Temperature:	4500 K
LED Lifespan:	50 000 Hrs
Colour Rendering Index:	>85
Switching:	Illuminated switch on lamp head.
Lamp body/Housing: (All plastic Parts)	Plastic moulded.
Arm:	Aluminium: Natural anodised.
Colour:	Lamp body - White, Arms - Natural anodised, cover caps offered in blue.
Mains Lead:	Moulded plug top, 1.5 m cable. (Wall) / Terminates at unit. (Ceiling)
Product Weight (Net):	2.0 Kg 9wall) / 7.5 Kg (Ceiling)
Dimensions (Packaging):	780 mm (H) x 200 mm (W) x 130 mm (D) *Folded (Wall) / 690 mm (H) x 260 mm (W) x 1012 mm (L) *Folded (Ceiling)
Operating Reach:	1200 mm (Wall) / 1540 mm (Ceiling)
Mounting Options:	Jack Russell-W (Wall Mounted) / Jack Russell-C (Ceiling Mounted)
Warranty:	2 Years
IP Rating:	54



OTHER CONFIGURATIONS:

