

ADVANCING HEALTHCARE TECHNOLOGY



PRODUCT CATEGORY:

MEDICAL AND SURGICAL LIGHTING

PRODUCT NAME:

LUVIS INTERNAL CAMERA - WIRELESS

PRODUCT DESCRIPTION:

OPERATING THEATRE LIGHT



Area of Use: Operating Theatres

Intended Use: Wireless transmission of the internal camera provides excellent video stream with a live video feed, enabling clinical teams to view the operating area from an external location.



LUVIS INTERNAL CAMERA - WIRELESS SPECIFICATIONS

CONTENT	SPECIFICATION		REMARK
	WIRE	WIRELESS	
Image Sensor:	1/2.8 type CMOS		
Zoom Ratio:	30 x Optical Zoom		
Image Point:	Approx. 2.13 Megapixels		
Min Object Distance:	10 mm - 1200 mm		
Video System:	1080p/60		
Video Output Signal:	3G HD-SDI SMPTE 424M & 292M		75Ω BNC
Wireless Frequency:	-	5.0 GHz	
Wireless Bandwidth:	-	20 GHz	
Wireless Mac Distance:	-	120 m	
Transmission Latency:	-	0.1 s	



OTHER PICTURES:

