

OUR TRUSTED LEGACY IS YOUR SUCCESS



PRODUCT CATEGORY:
VETERINARY LIGHTING

PRODUCT NAME:
CHIHUAHUA

PRODUCT DESCRIPTION:
HEADLAMP



Area of Use: Veterinary Rooms / Veterinary Theatres / Veterinary Clinics

Intended Use: The Chihuahua Headamp is designed to illuminate the surgical or examination site, improving visibility and aiding in accurate diagnosis and treatment. The device is intended to be used by Veterinarians / Veterinarian staff.

OUR TRUSTED LEGACY IS YOUR SUCCESS



CHIHUAHUA SPECIFICATIONS

SPECIFICATIONS:	
Rated Voltage:	AC100~240V 50/60Hz
Lamp & Electricity Usage:	LED 5W
Working Distance:	420 mm
Verification Distance:	200 mm
Colour Temperature:	4000 K - 5 500 K
Product Illumination:	≥ 30000 Lux
Charging Time:	≤ 5 h
Power Supply Time:	4h+4h
Weight:	330 g