

OUR TRUSTED LEGACY IS YOUR SUCCESS



PRODUCT CATEGORY:
MEDICAL LIGHTING

PRODUCT NAME:
HU-BM-HEADLAMP

PRODUCT DESCRIPTION:
HEADLAMP



Area of Use: Operating Theatre / Optometry / Dermatology

Intended Use: The Headamp is designed to illuminate the surgical or examination site, improving visibility and aiding in accurate diagnosis and treatment. The device is intended to be used by doctors.

OUR TRUSTED LEGACY IS YOUR SUCCESS



HU-BM-HEADLAMP SPECIFICATIONS

SPECIFICATIONS:	
Rated Voltage:	AC100~240V 50/60Hz
Lamp & Electricity Usage:	LED 5W
Working Distance:	420 mm
Verification Distance:	200 mm
Colour Temperature:	4000 K - 5 500 K
Product Illumination:	≥30000 Lux
Charging Time:	≤5h
Power Supply Time:	4h+4h
Weight:	330 g