

OUR TRUSTED LEGACY IS YOUR SUCCESS



PRODUCT CATEGORY:
VETERINARY LIGHTING

PRODUCT NAME:
DOBERMANN

PRODUCT DESCRIPTION:
VETERINARY PROCEDURE LIGHT



Area of Use: Veterinary Rooms / Veterinary Theatres / Veterinary Clinics

Intended Use: The Dobermann Veterinary Procedure Light is designed to provide illumination of the body of a animal for localized diagnosis and treatment. The device is intended to be used by Veterinarians / Veterinarian staff.

OUR TRUSTED LEGACY IS YOUR SUCCESS



DOBERMANN SPECIFICATIONS

SPECIFICATIONS:	
Multicoloured Light:	Yes
Central illuminance at 1000 mm:	100 000 Lux @ 1000 mm
Intensity Control from/to:	40 000 - 100 000 Lux - 5 steps
Adjustable Focus:	180 - 300 mm
Colour Temperature:	3800 K / 4 300 K / 4800 K - 3 steps
Reach:	1650 mm (From mount to centre of light source.)
Colour rendering index (Ra):	95
Number of LED's:	18
LED Lifespan:	>50 000 Hrs
Head size:	350 mm x 340 mm
Mounting Options:	Wall (Dobermann-W) /Ceiling Mounted (Dobermann-C)
Detachable Centre Grip:	Yes (x2)
Certification:	CE / FDA / CSA / UL
Additional Options:	
Wall Control Option:	N/A
Battery Backup Option:	N/A
Weight:	Single: 25 Kg



OTHER CONFIGURATIONS:



Wall Version

Images are used for illustration purposes only.