

ADVANCING HEALTHCARE TECHNOLOGY



PRODUCT CATEGORY:

MEDICAL AND SURGICAL LIGHTING

PRODUCT NAME:

LUVIS S200

PRODUCT DESCRIPTION:

MEDICAL PROCEDURE LIGHT



Product Codes: LUVIS S200

Area of Use: Normal Care to High Care Wards, Medical Consultation Rooms, Intensive Care Units, Trauma Units and Clinics.

Intended Use: The LUVIS S200 Medical Procedure Light is designed to provide illumination of the body of a patient for localized diagnosis and treatment. The light is not designed for use during surgery. It is designed for continuous operation and it must be possible to stop the treatment or diagnosis without danger to the patient if the light fails during use. The device is intended to be used by doctors / medical practitioners and nursing staff.



LUVIS S200 SPECIFICATIONS

SPECIFICATIONS:	
Multicoloured Light:	Yes
Central illuminance at 1000 mm:	100 000 Lux @ 1000 mm
Intensity Control from/to:	40 000 - 100 000 Lux - 5 steps
Adjustable Focus:	180 - 300 mm
Colour Temperature:	3800 K / 4 300 K / 4800 K - 3 steps
Colour rendering index (Ra):	95
Number of LED's:	18
LED Lifespan:	>50 000 Hrs
Head size:	350 mm x 340 mm
Mounting Options:	Wall (Single) /Mobile (Single) /Ceiling Mounted (Single or Dual)
Detachable Centre Grip:	Yes (x2)
Certification:	CE / FDA / CSA / UL
Additional Options:	
Wall Control Option:	N/A
Battery Backup Option:	N/A
Weight:	Single: 25 Kg



OTHER CONFIGURATIONS:



Images are used for illustration purposes only.